
Changing the Paradigm in Petroleomics with Comprehensive Two-Dimensional Gas Chromatography High Resolution Mass Spectrometry


Jonathan D. Byer¹, Sophia Schreckenbach², Gregg Tomy², Joseph E. Binkley¹, David E. Alonso¹, Lorne E. Fell¹ | ¹LECO Corporation, 3000 Lakeview Avenue, Saint Joseph, MI, United States I ²Department of Chemistry, University of Manitoba, Winnipeg, MB, Canada

Introduction


Petroleomics is the characterization of the chemical constituents of naturallyoccurring petroleum and crude oil using high resolution mass spectrometry.^{[1][2][3]} In addition to mass determination, petroleomic analysis sorts the chemical compounds into heteroatom class (nitrogen, oxygen, and sulfur), type (degree of unsaturation), and carbon number.^[4]


High resolution or ultra-high resolution MS is needed to resolve isobaric masses






GC×GC provides superior chromatographic separation for complex matrices






Is there any benefit in combining GC×GC and HRMS for Petroleomics? Mass Classification Chromatographic classification



Analytical Conditions

Primary Column: 60 m x 0.25 mm ID x 0.25 µm df (Rxi-17SilMS) **Sec. Column**: 0.6 m x 0.25 mm ID x 0.25 µm df (Rxi-5MS) **Sample:** 1 µL cold splitless **Column Flow**: 1.4 mL/min **Injection Temp**: 40° C – (720°C/min) – 300°C **Oven Temp**: 50°C (1.5 min) – (5°C/min) – 350°C (5.5 min); 67 min **Sec Oven:** +5°C **Modulator:** +15°C **Modulation Period:** 1.5 s MS Transfer Line Temp: 330°C Mass Range (EI): m/z 15–535 **Mass Range (CI):** m/z 45–535 (CH₄); m/z 60–535 (NH₃) Acquisition Rate: 150 spectra/s Resolution: 25,000 Source Temp: 250°C (EI); 230°C (CI)

MORE ANALYTES IDENTIFIED



- Compared to Low Resolution
- Increased mass resolution
- Chemical formula
- computation Unknown–unknown
- identification

Benefits of GC-HRT 4D

- Combination of GC×GC and accurate mass MS
- Highest acquisition rate, highest resolution, best mass accuracy of any GC×GC capable MS on the market
- Software for GC×GC-HRMS
- Fully integrated




This technique can be used to distinguish the C₃/SH₄ Mass Split

Theoretical mass split $C_3 vs SH_4 = 3.4 mDa$ ←→ C₁₈H₃r C18H14O 0.025 Da 3.5mDa **→** 246.1435 ${}_{17}^{\bullet}H_{10}S$ 246.0497 246.075 246.125 246.175 246.225 \leftrightarrow 1.0 Da ←→ 50 Da 200 100 150 200 250 300 350 400 450 500 M/Z

The integrated mass spectrum, based on spectrally deconvoluted EI data, of a dilbit sample for all of the identified species in the Kendrick plot above.


Three isobaric compounds with resolution requirements > 70 000 FWHH were distinguished by GC×GC-HRTOFMS operated at RP = 25 000 FWHH







Results

The XIC of m/z 246.14 shows three clusters of peaks. The accurate mass data for each of the peak clusters indicates there are three compound classes represented. The bottom two contour plots show an example of the $C_3 \SH_4$ mass split (0.0034 Da), an important feature in petroleomics. The figure also demonstrates isomer separations of each of the compound classes, a benefit not achievable with direct-infusion FT-MS.







Comprehensive two-dimensional GC coupled to HRTOFMS is a powerful tool for chemical characterization. This technique is ideal for petroleomics because it provides mass determination, heteroatom class, type, and carbon number information similar to direct high resolution mass spectrometry, and additionally gives well-defined chromatographic regions that facilitate structural elucidation and speciation.

[1] Marshall, Alan G.; Rodgers, Ryan P. (2004). Accounts of Chemical Research **37** (1): 53–59. [2] Marshall, A. G.; Rodgers, R. P. (2008). Proceedings of the National Academy of Sciences 105 (47). 18090-18095. [3] Ho, Yunju; Ahmed, Arif; Islam, Annana; Kim, Sunghwan (2014). Mass Spectrometry Reviews 34: 248 -263.

References

[4] Oliver C. Mullins; Eric Y. Sheu; Ahmed Hammami; Alan G. Marshall (8 November 2007). Asphaltenes, Heavy Oils, and Petroleomics. Springer. ISBN 978-0-387-68903-6.