
Analysis of Petroleum Products using Comprehensive Two-Dimensional Gas Chromatography (GC×GC) with Both Time of Flight MS and Flame Ionization Detectors

Introduction


Comprehensive Two-Dimensional Gas Chromatography (GC×GC) has proven to be an extremely valuable analytical technique for the petroleum industry due to its ability to substantially increase the chromatographic peak capacity beyond that of traditional single-dimension gas chromatography. Pairing GC×GC with Time-of-Flight Mass Spectrometry (TOFMS) provides unsurpassed characterization capabilities due to the separation power of GC×GC and the ability of TOFMS to provide rich data to deconvolution algorithms, which help unravel the complexity of difficult petroleum matrices. Petrochemical labs also often utilize flame ionization detectors (FID) for GC×GC to provide quantitative results via area percent calculations.

Configurations

Several possible configurations have been tested to develop a workflow which would allow either subsequent or simultaneous acquisitions of GC×GC TOFMS and FID data on the same GC×GC instrument.



Configuration 1: Two inlets are each equipped with a primary column (30 m x 0.25 mm x 0.25 µm Rxi-5MS) connected to a secondary column (2 m x 0.25 mm x 0.25 µm) One column set routes to the FID, while the other column set routes to the MS.




Configuration 3: A single primary column (30 m x 0.25 mm x 0.25 µm Rxi-5MS) is connected to two identical secondary columns (2 m x 0.25 mm x 0.25 µm) using a low dead-volume Y-connector. An additional restrictive guard column (0.1 mm) routed to the mass spectrometer accounts for the difference between atmospheric and vacuum outlet detectors.

Methods


	Pegasu	JS [®] GC-HT/4D	
Column Configuration			
Primary Column	RXI-5MS 30 m x 0.25 mm x 0.25 μm		
Secondary Column	Rxi-17Sil MS 2 m x 0.25 mm x 0.25 µm		
Temperature Program—Grob	1.00 min at 40°C, ramp 10°C/min to 150°C, ramp 30°C/min to 230°C and hold 5.00 min		
Temperature Program—Diesel	1.00 min at 40°C, ramp 5°C/min to 280°C and hold 5.00 min		
MS Parameters			
Transferline Temp (°C)	300	Acquisition Rate (s/s)	200
Source Temp (°C)	200	Mass Range	35-500
FID Parameters			
Data Collection Rate (Hz)	200		
Temperature (ºC)	300		

Christina N. Kelly, Jonathan D. Byer, David E. Alonso, Lorne E. Fell, Joseph E. Binkley | LECO Corporation, St. Joseph, MI


Configuration 2: Two independent column sets as described in Configuration (left) were attached to a single inlet via a dual-holed ferrule. One column set routes to the FID, while the other column set routes to the MS.



Configuration 4: One primary column (30 m x 0.25 mm x 0.25 µm Rxi-5MS) connected to a secondary column (2 m x 0.25 mm x 0.25 µm) is routed to a controlled flow splitter (CFT), which uses different lengths and diameters of restrictive tubing to evenly split analyte between the two detectors.



MS-to-FID correlation of peaks



Diesel Sample Analyte Tables: The bar graphs above show the FID Area Percent values for Benzene,1,2,3-Trimethyl, Naphthalene, 2-Methylnaphthalene, 1-Methylnaphthalene, and Triacontane for GCXGC (2D) runs in all four configurations.



Diesel Sample


Diesel Sample Contour Plot: The contour plots above show data acquired in Configuration 4. Diesel Sample analytes are identified by GCxGC-MS, with the corresponding FID peaks shown below. Retention times are closely matched, facilitating

Diesel Sample Selected Analytes

In every configuration with simultaneous acquisition, one of the chief difficulties was trying to match the flows between analyte streams, because of the difference between the atmospheric pressure of the FID and vacuum pressure from the MS. With the temperature ramp throughout the run, this also resulted in uneven splitting of analytes in Configurations 2 and 3. However, use of a controlled flow splitter (CFT) in Configuration 4 with properly calculated restrictor lengths and diameters going to each detector allowed for retention times to match within one modulation period. GCXGC allowed for enough chromatographic resolution of certain peaks so that pattern recognition could be used to correlate the MS and FID peaks in the more complex diesel sample. Area percent was automatically calculated for each peak using the S1 signal for FID data and the Deconvoluted Total Ion Chromatogram (DTIC) signal for the MS data.

Of the configurations presented here, Configuration 4 would be the preferred method for analysis of petroleum products. Configuration 4 allowed for close alignment of peak retention times, greatly facilitating data interpretation. In addition, the integrity of the chromatographic resolution was maintained throughout the run, with the analyte flow only split when very close to each detector, allowing each individual peak from the FID trace to be properly integrated so that more accurate area percent values could be calculated, which led to slightly lower values than in Configurations 1-3. While there was some loss of sensitivity due to the additional analyte flow split, the ability to acquire FID and MS data simultaneously greatly reduced acquisition time, making analysis with two detectors both simpler and more efficient.





Spectral Similarity: The table above shows the excellent average similarity scores of analytes of interest across all configurations, which allow for confident identification of MS peaks

Discussion

Conclusions