High Resolving Power Assessment for High Resolution Time-of-Flight Mass Spectrometry

Chang Samuel Hsu¹, Jav Lu¹, Viatcheslav Artaev², David Alonso², Kevin McNitt² and Joe Binklev² 1. Florida State University, Tallahassee, FL 32310; 2. LECO Corporation, St. Joseph, MI 49085

Introduction

Mass spectrometry has played an indispensable role in the characterization of complex petroleum and other fossil fuel resources at molecular level, especially the determination of elemental composition of components. The coupling of chromatography enhances its capabilities of detailed compositional analysis. The frequent absence of a molecular ion or novel fragment ion among hydrocarbons makes separation desirable. Gas chromatography (GC) provides high chromatographic resolving power to separate individual, even isomeric, species to facilitate extraction of informative mass spectral information. LECO Corporation's recent development of high resolution time of flight mass spectrometers capable of up to 100,000 resolving power, sub ppm mass accuracy and acquisition rates up to 200 sspectra per second have opened the door for pairing high resolution chromatography methods (UHPLC, GC, and GCxGC) with high resolution mass analysis. The work contained in this poster will focus on the use of gas chromatography separations paired with high resolution time of flight mass spectrometry on the LECO Pegasus® GC-HRT.

The analysis of petroleum biomarkers, such as hopanes and steranes. can provide insights into the quality of the crude oil as well as its source. Much information can be gleaned about a crude oil sample from its hopane and sterane content. Information such as source, age, maturity, and alteration. This information can also be a useful tool for oil to oil correlation as well as oil to source rock correlation. The ability to sufficiently resolve (mass spectrometrically and chromatographically) the key masses typical of hopane and sterane compounds present in Macondo Crude oil will be evaluated.


The objective of the current study is to compare and evaluate the utility of accurate mass and resolving power in the selective extraction of gualitative information on hopane and sterane content of a representative oil sample. In this instance a sample of Macondo crude is used. In addition, GCxGC data are also shown as a complementary form of selectivity.



Pegasus® GC-HRT with Patented Folded Flight Path FFP™ Mass Analyzer

Experimental

A sample of crude oil from the Macondo well was prepared in pentane prior to analysis. This sample was analyzed in each of the three modes. Three different analysis modes are available on the Pegasus GC-HRT: Nominal Resolution (R = 1,000), High Resolution (R = 25,000), and Ultra High Resolution (R = 50,000). The mass resolution is achieved by a mass analyzer with variable flight path length. The figure below illustrates the flight paths representative of each analysis mode.




The following GC-HRT Conditions Were Utilized:


GC Parameters

Sample


GC: Column Type: njection:	Agilent 7890 Rxi-5Sil MS (30m x0.25mm ID x 0.25 mm <i>df</i>) 1µL, Split 20:1 at 300°C
Oven:	60 °C (1 min.) to 300 °C at 30 °C /min, to 350 °C at 10 °C /min
Carrier Gas:	He, Constant Flow (1.00 mL/min.)
MS Parameters	
Spectrometer: on Source:	LECO PEGASUS GC-HRT® El
Polarity:	Positive

opectioniciei. Le	
Ion Source:	El
Polarity:	Positive
Flight Path:	Nominal, HR and UHR, and UHR
Spectral Acquisition:	6 spectra / second
m/z Range:	60 – 510 (N and HR); 90 – 340 (UHR)
m/z Calibration:	PFTBA





in nominal mode (left), interfering masses are not resolved. In the center and right examples (HR and UHR modes, respectively), additional interfering masses have been resolved.



m/z 217 (sterane fragment) above is observed in 3 different operation modes. It is apparent that in nominal mode (left), interfering masses are not resolved. In the center and right examples (HR and UHR modes, respectively), additional interfering masses have been resolved

- defining the appropriate analyte signal. Nominal mass resolution does not distinguish the chromatographically unresolved nominal isobars. Coupling with chromatography is necessary for resolving components of
- different compound classes with identical formula, i.e., same accurate mass.
- Coupling of comprehensive two-dimensional gas chromatography (GCxGC) with HRT shows promise and the potential to resolve additional isomers and provide mass resolution to distinguish components not sufficiently resolved.

References:

1. Beynon, J. H.; Williams, A. E. Mass and Abundance Tables for Use in Mass Spectrometry: Elsevier: Amsterdam, 1986.

- 2. Nikolaev, E. N.; Boldin, I. A.; Jert, R.; Bavkut, G. J. Am. Soc. Mass Spectrom., 2011, 22, 1125-1133, 3. Hsu. C. S. Energy Fuels. 2012. 26. 1169-1177.
- Song C.; Hsu, C. S.; Mochida, I. (eds.) Chemistry of Diesel Fuels. Philadelphia: Taylor & Francis, 2000. 5. Hsu, C. S. (ed.) Analytical Advances for Hydrocarbon Research. New York: Kluwer Academic/Plenun ublishers, 2003 (ISBN 0-306-47476-X)