


What Else Is In My Dioxin Sample? High Performance Time-Of-Flight Mass Spectral Analysis Of Environmental Samples

Slava Artaev¹, Eric J. Reiner², Peter Gorst-Allman³, Joe Binkley¹, and David E Alonso¹

1. LECO Corporation, USA 2. University of Toronto 3. LECO Africa


Introduction

Polychlorinated dibenzo-*p*-dioxins and furans are some of the most toxic materials known to man. They are relatively inert, fat soluble halogenated organic compounds that bioaccumulate in the environment. These persistent organic pollutants (POPs) have been linked to adverse cancer, reproductive and endocrine effects. In addition to dioxins, polychlorinated biphenyls (PCBs) and polybrominated diphenyl ethers (PBDEs), environmental samples may contain polycyclic aromatic hydrocarbons (PAHs), halogenated PAHs, or numerous other industrial products and undiscovered hazardous compounds. These untargeted compounds may also be toxic and can be present in much higher concentrations than the target dioxins.



Persistent Organic Pollutants

High resolution TOFMS can be used for the comprehensive analysis of complex environmental samples. The advantages of high performance TOFMS for POP analysis are rapid acquisition of spectral data across a wide mass range with a minimum reduction in sensitivity, high resolving power to minimize interferences and high mass accuracy for robust elemental composition determinations. While other MS systems rely on selected ion monitoring to achieve the necessary detection limits for environmental studies, high resolution TOFMS can provide sensitivity while acquiring full mass range spectra allowing for identification of not only targeted compounds, but newly emerging environmental contaminants.



High Performance TOF MS (left) and its Folded Flight Mass Analyzer (Right)

Experimental

Standards and Samples:


Calibration standards (CS1 – CS5) were purchased from Wellington Laboratories (Guelph, Canada). Sediment and fly ash samples containing both PCDDs and PCDFs, as well as several classes of POPs were obtained from the EPA in Taiwan.

Instrument Parameters:


- GC: Agilent 7890
 - Column Type: Restek Rxi-5 (30 m, 0.25 mm ID, 0.25 μm d)
 - Inj. Temp.: 250 C
 - Injection: Splitless, 2 μL
 - Oven: 120 C (1)→220 C(20 C/min)→240 C(2.0 C/min)→250 C (1.0 C/min)→260 C(5.0 C/min)→265 C(1.0 C/min)
 - Carrier Gas: He, 1.0 mL/min constant flow
- MS: Pegasus GC-HRT
 - Source Temp.: 250 °C
 - EI: 40 eV
 - Mass Range: 160 - 510 m/z, High Resolution Mode

Representative Standard Data: CS5

Mass Spectra (OCDD, OCDF):




Relative Abundance Values (OCDD, OCDF)



Data CS3 Standard

Extracted Ion Chromatograms (XICs) -- Five Replicate Injections:



Mass Accuracy Values: Native PCDF/Ds

Name	Formula	Calculated M _r	Observed M _r	Mass Accuracy (ppm)
TCDF N	C12H4OC4	303.90108	303.90123	0.51
TCDD N	C12H2OC4	319.89599	319.89556	-1.36
PCDF N	C12H3OC5	337.86211	337.86202	-0.26
PCDF N	C12H3OC5	337.86211	337.86183	-0.81
PCDD N	C12H2OC5	353.85702	353.85716	0.41
HxCDF N	C12H2OC6	371.82313	371.82306	-0.20
HxCDF N	C12H2OC6	371.82313	371.82354	1.09
HxCDF N	C12H2OC6	371.82313	371.82302	-0.30
HxCDD N	C12H2OC6	387.81805	387.81810	0.14
HxCDD N	C12H2OC6	387.81805	387.81817	0.32
HxCDD N	C12H2OC6	387.81805	387.81830	0.65
HxCDF N	C12H2OC6	371.82313	371.82348	0.94
HxCDF N	C12H2OC7	405.78416	405.78394	-0.55
HxCDD N	C12H2OC7	421.77908	421.77908	0.00
HxCDF N	C12H2OC7	405.78416	405.78416	0.00
OCDD N	C12O2C8	455.74010	455.74044	0.74
OCDF N	C12OC8	439.74519	439.74450	-1.56

Ave. = 0.57 ppm


Mass Accuracy Values: Labeled PCDF/Ds

Labeled Compound	Calc M _r	Obs M _r	Mass Accuracy (ppm)	Calc M _r 2	Obs M _r 2	Mass Accuracy (ppm)
TCDF L	315.94134	315.94136	0.08	317.93839	317.93851	0.39
TCDD L	331.93625	331.93650	0.76	333.93330	333.93331	0.03

Labeled Compound	Calc M _r 2	Obs M _r 2	Mass Accuracy (ppm)	Calc M _r 4	Obs M _r 4	Mass Accuracy (ppm)
PCDF L	351.89941	351.89951	0.28	353.89646	353.89646	-0.01
PCDF L	351.89941	351.89943	0.05	353.89646	353.89653	0.19
PCDF L	367.89433	367.89474	1.12	369.89134	369.89154	0.44
HxCDF L	385.86044	385.86054	0.26	387.85749	387.85761	0.31
HxCDF L	385.86044	385.86034	-0.26	387.85749	387.85787	0.98
HxCDF L	385.86044	385.86057	0.34	387.85749	387.85732	-0.44
HxCDD L	401.85535	401.85543	0.19	403.85240	403.85215	-0.63
HxCDD L	401.85535	401.85558	0.56	403.85240	403.85259	0.46
HxCDD L	401.85535	401.85555	0.49	403.85240	403.85252	0.29
HxCDF L	385.86044	385.86034	-0.26	387.85749	387.85754	0.13
HxCDF L	419.82147	419.82138	-0.21	421.81852	421.81864	0.29
HxCDD L	435.81638	435.81667	0.66	437.81343	437.81322	-0.49
HxCDF L	419.82147	419.82162	0.36	421.81852	421.81878	0.62
OCDD L	469.77741	469.77746	0.11	471.77446	471.77413	-0.70


Ave. = 0.39 ppm

Calibration Results:



Results: Fly Ash


TIC and XIC for sample:




Target Analyte Finding -- 76 Dioxins:

Peak #	Name	R.T. (s)	Quant (μg/g)	Area	Height
1	TCDF N	607.5	231	45522	6029
2	TCDF N	617	146	27600	2969
3	TCDF N	626	70	9492	1897
4	TCDF N	636	1500	59957	36533
5	TCDD N	645	85	11796	1700
6	TCDD N	650.5	98	13372	1935
7	TCDF N	655	341	94652	6811
8	TCDF N	671	1193	50375	23865
9	TCDF N	681.5	419	87513	3875
10	TCDD N	684.5	245	93572	4901
11	TCDD N	700.5	71	15029	1425
12	TCDF N	701	866	144534	8660
14	TCDD L	711	585	152418	15095
15	TCDF N	714	1355	355656	26506
16	TCDD N	720.5	146	8007	2156
18	TCDF N	736	47	5978	1149
19	TCDD N	751.5	25	7465	506
20	TCDF N	761.5	100	30071	1993
21	PCDF N	800	394	181936	7392
22	PCDF N	824	681	104802	13618
23	PCDF N	841.5	64	1745	1378
24	PCDF N	846	2062	359600	21488
25	PCDF N	856	54	5874	1076
27	PCDF N	864.5	672	67054	8726
28	PCDF N	867	324	69060	6483
29	PCDF N	876.5	164	20765	3476
30	PCDD N	879.5	305	71588	6097
31	PCDF N	881.5	147	34460	2442
32	PCDF N	884.5	176	48001	3429
33	2378-PCDF N-2	909	2049	440715	20429

What else is in my sample?




Mass Spectra: PAHs




Current Work

Current work focuses on the use of a pulsed electron ionization source for sensitivity enhancement, and improving chromatographic resolution by using longer, thin film columns (e.g., Rtx-Dioxin-2, 40m x 0.18 mm x 0.18 μm).


XIC 2378-TCDF and 2378-TCDD (CSL – CS1)



Chromatographic Resolution: CS3WT Verification Standard



XIC sediment sample – 110 Dioxins:



Conclusion

High performance TOFMS provides a comprehensive profile of samples in a single acquisition. Its high resolving power and excellent mass accuracy values make it a practical choice for the analysis of complex environmental samples.