

Non-targeted analysis of contaminants in dust using multidimensional chromatography, high-resolution mass spectrometry and advanced software tools



David E. Alonso and Joseph E. Binkley; LECO Corporation, Saint Joseph, Michigan, USA


Introduction

- Dust is complex: Numerous, chemically diverse constituents in a wide concentration range
- Dust can be used to estimate human exposure to contaminants (e.g., Pesticides, PAHs, flame retardants)
- Analysis of dust has primarily been conducted using targeted methods and a variety of instrumental techniques

Objectives

- Untargeted characterization of dust extracts using GCxGC and a novel multimode source with high resolution time-of-flight mass spectrometry
- To perform quantitative analysis of Polychlorinated Paraffins (PCPs) in dust

Analytical Platform



Mass Accuracy: 1 PPM
Acquisition Speed: 200 sps
Enhanced Chromatographic Resolution: GCxGC
Resolution: Up to 50,000
Multimode Source (MMS): EI, PCI & ECNI

Figure 1. HRT+ 4D and Multimode Source.

Dust Sample Preparation

- Dust: 1) NIST SRM 2585, 2) Office Sample, and 3) Household Sample
- Extraction:




Figure 2. Dust general extraction procedure.


Instrument Acquisition Parameters

Gas Chromatograph	Agilent 7890B with LECO Dual Stage Quad Jet Modulator
Injection	2µL liquid injection, Splitless, 70°C to 300°C at 500°C/min
Carrier Gas	He @ 1.0 mL/min, Corrected Constant Flow
Primary Column	HP-5MS UI, 30 m x 0.25 mm i.d. x 0.25 µm
Secondary Column	BPX-50, 0.60 m x 0.10 mm x 0.10 µm
Temperature Program	80°C (1 min) ramp 20°C/min to 140°C (Hold 30 mins), then ramp to 200°C at 10°C/min, and ramp 1.5°C/min to 300°C (Hold 5 mins). Secondary oven maintained +5 °C relative to primary oven
Modulation Period	8.0 seconds; modulator maintained +15 °C relative to secondary oven
Transfer Line	300 °C
Mass Spectrometer	LECO Pegasus® HRT+ 4D
Source Temperature	EI, 250°C; PCI, ECNI, 165°C
Acquisition Mode	High Resolution, R ≥ 25,000 for m/z 219, Mass Accuracy ≤ 1 ppm
Ionization	EI, PCI/ECNI (Reagent Gas = CH ₄)
Mass Range (m/z)	EI 50-1000; PCI 60-1000; ECNI 30-1000
Acquisition Rate	MMS 125 sps

Table 1. HRT+ 4D instrument acquisition parameters.

- Hydrocarbons
- Acids
- Aromatics
- Amines
- Alcohols
- Aldehydes
- Ketones
- Phenols
- Fatty Acids
- Sterols
- Phosphates
- More...

NIST SRM: EI, ECNI & PCI-HRT+ 4D Data



- Hard ionization (Fragmentation = More Information)
- Non-selective (Universal)
- Good sensitivity
- Soft ionization (Less Fragmentation)
- Selective (Electron Rich Compounds)
- Increased sensitivity

- Soft ionization (Less Fragmentation)
- Selective (Electron Rich Compounds)

Figure 3. EI, ECNI, and PCI contour plot displaying NIST SRM 2585 data collected using the Multimode Source and GCxGC-TOFMS.

Advantages of GCxGC-HRT+ 4D Data




Figure 4. EI plot displaying some of the major compounds in SRM 2585.

NIST SRM 2585: Representative Compounds

Name	Formula	R.T. (s)	Similarity	Name	Formula	R.T. (s)	Similarity	Name	Formula	R.T. (s)	Similarity	
Benzyl chloride	C ₇ H ₇ Cl	488.1952	873	0.87	2,5-Dichloroaniline	C ₆ H ₄ Cl ₂ N	1872.5501	750	2,5-Dichloroaniline	C ₆ H ₄ Cl ₂ N	1872.5501	750
Benzyl bromide	C ₇ H ₇ Br	560.2263	866	0.84	2-Chlorophenyl isocyanate	C ₇ H ₅ ClN	1704.4472	874	2-Chlorophenyl isocyanate	C ₇ H ₅ ClN	1704.4472	874
3-Chlorophenol	C ₆ H ₄ ClO	664.2784	896	1.32	4-Chloroaniline	C ₆ H ₄ ClN	696.3300	933	4-Chloroaniline	C ₆ H ₄ ClN	696.3300	933
4-Chloroaniline	C ₆ H ₄ ClN	704.3184	836	1.13	1-Chlorobromobenzene	C ₆ H ₄ BrCl	1554.3556	949	1-Chlorobromobenzene	C ₆ H ₄ BrCl	1554.3556	949
1,2-Dichlorobenzene	C ₆ H ₄ Cl ₂	2096.3768	873	N/A	2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	2488.2976	783	2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	2488.2976	783
2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	2778.3064	836	N/A	2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	2778.3064	836	2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	2778.3064	836
2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	2904.4100	824	1.22	2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	2904.4100	824	2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	2904.4100	824
2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	4384.6528	833	0.34	2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	4384.6528	833	2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	4384.6528	833
2,2',4,4',5,5'-Pentabromodiphenyl ether	C ₁₂ H ₆ Br ₅ O	5056.2296	892	-0.45	2,2',4,4',5,5'-Pentabromodiphenyl ether	C ₁₂ H ₆ Br ₅ O	5056.2296	892	2,2',4,4',5,5'-Pentabromodiphenyl ether	C ₁₂ H ₆ Br ₅ O	5056.2296	892

Ave. 880/1000

Table 2. Representative compounds in SRM 2585.

NIST SRM 2585: Halogenated Organic Compounds




Figure 5. EI plot displaying halogenated compounds in SRM 2585.

Name	Formula	R.T. (s)	Similarity	Mass Accuracy (ppm)
Benzyl chloride	C ₇ H ₇ Cl	488.1952	873	0.87
Benzyl bromide	C ₇ H ₇ Br	560.2263	866	0.84
2-Chlorophenyl isocyanate	C ₇ H ₅ ClN	696.3300	933	0.31
3-Chlorophenol	C ₆ H ₄ ClO	664.2784	896	1.32
4-Chloroaniline	C ₆ H ₄ ClN	704.3184	836	1.13
p-Tolylbenzyl chloride	C ₉ H ₉ Cl	798.2604	884	0.10
p-Tolylbenzyl bromide	C ₉ H ₉ Br	792.2704	933	0.21
o-Bromostyrene	C ₈ H ₇ Br	912.3528	862	-0.23
3,4-Dichlorophenyl isocyanate	C ₇ H ₃ Cl ₂ N	944.3358	820	1.93
2,5-Dichloroaniline	C ₆ H ₄ Cl ₂ N	1352.7124	933	-0.25
1-Chlorobromobenzene	C ₆ H ₄ BrCl	1554.3556	949	-0.80
Dodecane, 1-bromo-	C ₁₂ H ₂₅ Br	2096.3768	873	N/A
2,3,4-Trichloroacetophenone	C ₈ H ₅ Cl ₃ O	2488.2976	783	N/A
TCPP	C ₁₀ H ₇ Cl ₃ O ₃ P	2632.3302	893	N/A
TCBP	C ₁₀ H ₇ Cl ₃ O ₂ P	2704.3008	938	N/A
TCPP	C ₁₀ H ₇ Cl ₃ O ₃ P	2778.3064	836	N/A
Chlorophene	C ₁₀ H ₇ ClO	2904.4100	824	1.22
TCBP	C ₁₀ H ₇ Cl ₃ O ₂ P	3912.4504	902	N/A
2,3,4,5-Tetrabromodiphenyl ether	C ₁₂ H ₆ Br ₄ O	4384.6528	833	0.34
2,2',4,4',5,5'-Pentabromodiphenyl ether	C ₁₂ H ₆ Br ₅ O	5056.2296	892	-0.45

Table 3. SRM 2585 selected halogenated compounds.




Figure 6. Peak true and library EI mass spectra for TDCPP and BDE-99 in SRM 2585.

ECNI: Locating Trace POPs in NIST, Office & House Dust




Figure 7. Target Analyte Finding (TAF) method for POPs in dust samples.

Figure 8. ECNI plot displaying POPs in SRM 2585.

Target POP	R.T. (s)	NIST (area)	Office (area)	House (area)
C ₁₂ Cl ₆	2536, 2.872	66777	10364	485688
CN Cl ₄	2848, 3.304	13710		
CN Cl ₂	2920, 3.568	143301		
Chlordane Isomer 1	3136, 3.552	62129		
Heptachlor 1	3232, 4.128	15887		
CB Cl	3272, 4.192	21645	8827	
Heptachlor 2	3336, 4.360	215913		
trans-Chlordane	3360, 4.112	1730979	3025	11470
CB Cl ₂	3388, 4.552	10321	1454	
Chlordane Isomer 2	3392, 4.096	265412		
CB SCI	3400, 4.224	54160	50338	
Nonachlor Isomer 1	3432, 3.848	9577		
cis-Chlordane	3440, 4.176	624489		
Nonachlor	3464, 3.888	596640		3023
CB SCI-2	3584, 4.720	44766	55571	
CB SCI	3640, 4.536	23649	3324	
CB Cl ₂	3648, 4.520	23937		
CB SCI-3	3712, 4.712	304119	305058	
Nonachlor Isomer 3	3784, 4.656	98500		
CB Cl ₃	3800, 4.576	44400	46349	

Table 4. TAF Processing results for POPs in SRM 2585, office, and household dust.

Polychlorinated Paraffins: Cloud Quant




Figure 9. ECNI plot for PCPs in a calibration standard mix. Each group (cloud) region represents PCPs for which the number of carbon + chlorine atoms is constant.




Figure 10. ECNI Plot for Group 4 PCPs with formulas C₁₁H₁₇Cl₇ and C₁₂H₂₀Cl₆. Mass accuracy values for the most abundant isotope in each cluster was -1.14 & 0.08 ppm respectively.




Figure 11. ECNI Plot for SCCPs and MCCPs in the NIST SRM and a household dust sample.

Summary

- HRT & MMS technology are valuable tools for the analysis of complex samples
- The highly ordered, comprehensive contour plots can be used to target trace POPs in dust with pinpoint accuracy
- POPs in Dust: SRM 2585 >> Office Dust > Household Samples
- Quantitative analysis of PCPs in dust samples was conducted using GCxGC-ECNI-TOFMS