Jonathan D. Byer<sup>1</sup>, Grazina Pacepavicius<sup>2</sup>, Peter V. Hodson<sup>3</sup>, Claude Belpaire<sup>4</sup>, David Alonso<sup>1</sup>, Joe Binkley<sup>1</sup>, Mehran Alaee<sup>2</sup> <sup>1</sup>LECO Corporation, Saint Joseph, MI, USA | <sup>2</sup>Aquatic Contaminants Research Division, Environment Canada | <sup>4</sup>Research Institute for Nature and Forest, Groenendaal-Hoeilaart, Belgium

## Introduction

Time-of-flight mass spectrometers (TOFMS) have gained popularity over scanning instruments for non-target and post-target analysis because:

- Full mass range spectra are acquired
- High acquisition rates can be achieved (>200 spectra/s)
- Minimal mass bias
- Spectral deconvolution using modern software

In this study, we used non-target analysis to determine the occurrence of lipophilic halogenated contaminants in American and European eel (Anguillid sp.) not previously identified by traditional target analysis.

**Methods** 

- Endangered Wildlife in Canada.
- Species



# that are not routinely monitored.



- pollutants
- seconds.



Mass spectra collected in El mode with mass range m/z = 45-1000 and resolving power >25, 000 for m/z = 218.9856.

## Background

• American eel are listed as a threatened species by the Committee on the Status of • European eel are considered critically endangered by the IUCN Red List of Threatened

Dekker, W., *Fisheries* **2003**, 28, (12), 28–30.

• Bioaccumulation of halogenated chemical contaminants has lead to a reduction in spawner quality and may be related to recruitment decline.

Objective: To identify halogenated contaminants in eels using high resolution TOFMS

• Fish from these sample locations are historically contaminated with polychlorinated biphenyls (PCBs; RT 1000-1800s), and other legacy halogenated persistent organic

• The AICs show that there are spatial differences in the chemical contamination profile of eels among these sampling sites, particularly between retention times 1200-1400



contaminants and their breakdown products.



# Non-Target and Post-Target Analysis of Emerging Halogenated Contaminants in American and European Eels

(379.7449+343.7868, respectively) for all seven sampling sites.

Name	Formula	Expected	Observed	Mass Delta (Da)	Mass Accuracy (ppm)
Pentabromobenzene	C <sub>6</sub> HBr <sub>5</sub>	471.5949	471.5948	-0.0001	-0.27
Octachlorostyrene	C <sub>8</sub> Cl <sub>8</sub>	379.7444	379.7441	-0.0003	-0.80
Heptachlorostyrene	C <sub>8</sub> HCl <sub>7</sub>	343.7863	343.7859	-0.0004	-1.18
Pentachloroanisole	$C_7H_3CI_5O$	279.8592	279.8591	-0.0001	-0.24
Trichlorophenol	C <sub>6</sub> H <sub>3</sub> Cl <sub>3</sub> O	195.9244	195.9245	-0.0001	-0.08