

Comparative PCDD/F Analysis with GC-HRMS, GC-HRTOFMS and GCxGC-TOFMS Discovery of Compounds Not Found in Environmental Analysis Guided by EPA 1613B

Peter Gorst-Allman¹, David Alonso², Jayne de Vos³, Jack Cochran⁴ and Eric Reiner⁵

¹LECO Africa, Kempton Park, Gauteng, RSA; ²LECO Corporation, Saint Joseph, MI, USA; ³National Metrology Institute of South Africa, Pretoria, Gauteng, RSA, ⁴Restek Corporation, Bellefonte, PA, USA, and ⁵Ontario Ministry of the Environment, Toronto, Canada

Introduction

Globally, chlorinated dioxins and furans are analysed by expensive gas chromatography – high resolution mass spectrometry (GC-HRMS) using EPA Methods such as 8290A or 1613B, or local methodology based on these. To achieve the sensitivity required by regulatory methods, the HRMS approach uses selected ion monitoring (SIM). This is a target molecule technique, and requires a list of the selected analytes before running the analysis to set up the experimental parameters. Only the targeted compounds are located – all other information about the sample is forfeited. Using high res systems to do SIM is also expensive and requires a high degree of operator sophistication.

Time of Flight (TOF) mass spectrometers do not use SIM. All analyses provide full range mass spectra for all the analytes in the sample.


This is an important advantage, and provided the mandated sensitivities can be obtained using TOF, we should then have a method which can locate and identify all components of the sample, while also being able to quantify targeted PCDDs and PCDFs, all in one run.

This applies both to the Pegasus 4D (GCxGC-TOFMS) and to the Pegasus HRT (High Resolution TOFMS).

Determine if GCXGC-TOFMS and GC-HRT has the sensitivity to calibrate 2378-TCDD and 2378-PCDF down to 0.5 pg/μl as required by EPA Method 1613, and to analyse samples at this level.

Run the relevant ¹³C labelled standard set for 1613B
Use these results to obtain response factors for the PCDD/Fs
Run the sample set and calculate quantitative values for the PCDD/Fs
Data mine the samples for other POPs which may be present in the samples

Methods



Analysis Conditions GC-HRMS

Inlet 280°C
Constant Flow He at 0.8 ml/min
40m x 0.18mm x 0.18μm DB-5
140°C (1min), 52°C/min to 200°C, 2.9°C/min to 235°C (3 min), 3°C/min to 267°C (3 min), 7°C/min to 310°C (Hold until OCDD elutes)

Transfer Line 280°C
HRMS, EI+, SIM, 35 eV

Analysis Conditions GC-HRT

Inlet 280°C
Corrected Constant Flow He at 1 mL/min
40m x 0.18mm x 0.18μm Rtx-Dioxin2 Arylene backbone-modified siloxane column selective for coplanars
140°C (1min), 50°C/min to 200°C, 3°C/min to 260°C, 1°C/min to 280°C, 6°C/min to 310°C (5 min)
Transfer Line 300°C
HRT, EI+, 140 to 520 u, 3 spectra/s, 50 eV

Analysis Conditions GCxGC-TOFMS

40m x 0.18mm x 0.18μm Rtx-Dioxin2 Arylene backbone-modified siloxane column selective for coplanars
120°C (2min), 20°C/min to 200°C, 5°C/min to 320°C (3 min)
1.0m x 0.15mm x 0.15μm Rx-17Sil MS 50% phenyl (silylarylene) / 50% methyl type siloxane
+ 5°C offset from primary column
Corrected Constant Flow He at 1.4 mL/min (Speed optimized flow)
Thermal modulation, 2.0 sec (Hot time 0.70 sec, Cool time 0.30 sec)
TOFMS, EI+, 45 to 750 u, 100 spectra/s

Quantitative Results

| Sample 1 | 2378-TCDF | 2378-TCDD | 12378-PCDF | 23478-PCDF | 12378-PCDD | 123678-HxCDF | 123478-HxCDF | 123789-HxCDF | 123478-HxCDD |
|------------------|--------------|--------------|--------------|---------------|---------------|--------------|--------------|--------------|--------------|
| GC-HRT | 3.5 | 3.0 | 2.7 | 4.5 | 1.46 | 2.2 | 1.7 | ND | ND |
| GC-HRMS | 3.7 | 2.6 | 1.7 | 5.6 | 1.1 | 5.4 | 1.7 | 1.5 | 0.81 |
| GCxGC-TOFMS (x5) | ND | ND | ND | ND | ND | 11 | ND | ND | ND |
| | 3.5 | 5.6 | 6.6 | 12 | ND | 9.8 | 5.7 | 6.5 | ND |
| | 123678-HxCDD | 123789-HxCDD | 234678-HxCDF | 1234678-HxCDF | 1234789-HxCDF | OCDF | OCDD | | |
| GC-HRT | 3.1 | ND | ND | 12 | 33 | ND | 19 | 150 | |
| GC-HRMS | 1.7 | 1.6 | 0.60 | 15 | 24 | 1.5 | 24 | 170 | |
| GCxGC-TOFMS (x5) | ND | ND | ND | 31 | 16 | ND | 160 | | |
| | 6.3 | ND | ND | 31 | 16 | 31 | 170 | | |

| Sample 2 | 2378-TCDF | 2378-TCDD | 12378-PCDF | 23478-PCDF | 12378-PCDD | 123678-HxCDF | 123478-HxCDF | 123789-HxCDF | 123478-HxCDD |
|-------------|--------------|--------------|--------------|---------------|---------------|--------------|--------------|--------------|--------------|
| GC-HRT | 20 | 28 | 21 | 21 | 2.8 | 76 | 32 | 23 | 10 |
| GC-HRMS | 26 | 32 | 18 | 19 | 2.8 | 85 | 35 | 17 | 3.7 |
| GCxGC-TOFMS | 22 | 48 | 18 | 20 | ND | 90 | 40 | 32 | ND |
| | 123678-HxCDD | 123789-HxCDD | 234678-HxCDF | 1234678-HxCDF | 1234789-HxCDF | OCDF | OCDD | | |
| GC-HRT | 13 | 14 | 7.9 | 280 | 170 | 55 | 2400 | 1200 | |
| GC-HRMS | 9.3 | 6.0 | 2.3 | 300 | 160 | 47 | 2600 | 1900 | |
| GCxGC-TOFMS | ND | ND | ND | 410 | 230 | 120 | 2400 | 1500 | |

| Sample 3 | 2378-TCDF | 2378-TCDD | 12378-PCDF | 23478-PCDF | 12378-PCDD | 123678-HxCDF | 123478-HxCDF | 123789-HxCDF | 123478-HxCDD |
|-------------|--------------|--------------|--------------|---------------|---------------|--------------|--------------|--------------|--------------|
| GC-HRT | 48 | 55 | 22 | 140 | 10 | 220 | 37 | 42 | 17 |
| GC-HRMS | 46 | 69 | 14 | 180 | 3.1 | 240 | 21 | 14 | 5.2 |
| GCxGC-TOFMS | 40 | 73 | 32 | 150 | 21 | 210 | 46 | 41 | 18 |
| | 123678-HxCDD | 123789-HxCDD | 234678-HxCDF | 1234678-HxCDF | 1234789-HxCDF | OCDF | OCDD | | |
| GC-HRT | 18 | 28 | 13 | 620 | 240 | 51 | 980 | 790 | |
| GC-HRMS | 14 | 8.6 | 1.5 | 990 | 230 | 17 | 1100 | 1300 | |
| GCxGC-TOFMS | 41 | 32 | 36 | 950 | 210 | 58 | 980 | 990 | |

| Sample 4 | 2378-TCDF | 2378-TCDD | 12378-PCDF | 23478-PCDF | 12378-PCDD | 123678-HxCDF | 123478-HxCDF | 123789-HxCDF | 123478-HxCDD |
|-------------|--------------|--------------|--------------|---------------|---------------|--------------|--------------|--------------|--------------|
| GC-HRT | 52 | 3.2 | 120 | 39 | 5.1 | 280 | 190 | 8.9 | 5.1 |
| GC-HRMS | 59 | 3.9 | 140 | 55 | 4.7 | 330 | 210 | 28 | 4.1 |
| GCxGC-TOFMS | 48 | 5.8 | 130 | 36 | 14 | 210 | 210 | 75 | 23 |
| | 123678-HxCDD | 123789-HxCDD | 234678-HxCDF | 1234678-HxCDF | 1234789-HxCDF | OCDF | OCDD | | |
| GC-HRT | 13 | 20 | 86 | 710 | 65 | 380 | 6500 | 170 | |
| GC-HRMS | 9.9 | 8.1 | 50 | 1000 | 70 | 470 | 5200 | 220 | |
| GCxGC-TOFMS | 14 | 18 | 100 | 1000 | 66 | 460 | 5100 | 190 | |

| Sample 5 | 2378-TCDF | 2378-TCDD | 12378-PCDF | 23478-PCDF | 12378-PCDD | 123678-HxCDF | 123478-HxCDF | 123789-HxCDF | 123478-HxCDD |
|-------------|--------------|--------------|--------------|---------------|---------------|--------------|--------------|--------------|--------------|
| GC-HRT | 40 | ND | 35 | 24 | ND | 97 | 35 | 23 | ND |
| GC-HRMS | 50 | 0.65 | 40 | 31 | 1.9 | 120 | 48 | 20 | 2.9 |
| GCxGC-TOFMS | 43 | ND | 43 | 30 | ND | 120 | 47 | 29 | ND |
| | 123678-HxCDD | 123789-HxCDD | 234678-HxCDF | 1234678-HxCDF | 1234789-HxCDF | OCDF | OCDD | | |
| GC-HRT | ND | 14 | 17 | 290 | 71 | 72 | 1100 | 370 | |
| GC-HRMS | 5.8 | 3.9 | 4.7 | 420 | 75 | 79 | 1200 | 550 | |
| GCxGC-TOFMS | 10 | 15 | 38 | 400 | 78 | 90 | 1200 | 480 | |

| Sample 6 | 2378-TCDF | 2378-TCDD | 12378-PCDF | 23478-PCDF | 12378-PCDD | 123678-HxCDF |
| --- | --- | --- | --- | --- | --- | --- |