

# Fast Pesticide Residue Analysis Using a Novel Benchtop Time-of-Flight Mass Spectrometer

Jonathan D. Byer, Todd Richards, Joe Binkley, and Lorne Fell | LECO Corporation, Saint Joseph, MI USA

## Introduction

- Primary objective was to maximize throughput without introducing a significant number of additional chromatographic coelutions not separable by deconvolution, while transitioning from a standard 30 m column to a 15 m column.
- LECO's proprietary NonTarget Deconvolution™ (NTD)™ was leveraged to maintain exceptional peak fidelity of the qualitative analysis, while the Target Analyte Find data processing feature was utilized for robust quantitation and to establish method detection limits for organonitrogen pesticides spiked in QuEChERS strawberry extracts.


## GC Method for Fast Run Time

A bulk extract of strawberries purchased from a local grocery store was generated using methods described elsewhere—  
<http://www.restek.com/pdfs/GNAN1097A.pdf>.



A dilution series from 5000 ng/g to 0.10 ng/g of GC Multiresidue Pesticide Mix #5 (Restek) in the bulk extract was prepared in duplicate for GC-MS analysis, as well as a raw extract unfortified and without cleanup to investigate the occurrence of incurred pesticides. The instrument conditions used are shown in the table below.


|                               |  |
|-------------------------------|--|
| <b>Gas Chromatograph</b>      | <b>LECO L-PAL3 Autosampler, Agilent 7890B GC</b>                       |
| <b>Injection</b>              | 1 µL pulsed splitless, 20 psi for 0.7min, GC injector @ 250°C          |
| <b>Carrier Gas</b>            | He @ 2.0 ml/min, Constant Flow   |
| <b>Column</b>                 | Rxi-5ms, 15 m x 0.25 mmID x 0.25 µm df (Restek, Bellefonte, PA, USA)   |
| <b>Oven Program</b>           | 70 °C (0.7 min), to 150 °C @ 60 °C/min, to 330°C @ 30 °C/min (0.5 min) |
| <b>Transfer Line</b>          | 300 °C   |
| <b>Mass Spectrometer</b>      | <b>LECO Pegasus® BT</b>  |
| <b>Ion Source Temperature</b> | 225 °C   |
| <b>Mass Range</b>             | 45-650 m/z   |
| <b>Acquisition Rate</b>       | 10 spectra/s   |



Analytical ion chromatograms (AIC) of a fortified strawberry QuEChERS extract; (Red, top) for analytes reported using NonTarget Deconvolution (NTD) and (Blue, bottom) analytes reported using Target Analyte Find. The intense peaks in the Blue trace represent analytes that were incurred in the matrix in addition to the 10 ppb spiked. Spiked pesticides and their LODs are shown in Table I.

## Excellent Linearity and Low LODs in Matrix

Calibration curves for four representative pesticides in strawberry.




**Table I. Calibration curve linearity and LODs for organonitrogen pesticides in strawberry from 0.10 to 5000 ng/g using Target Analyte Find on a Rxi-5ms, 15 m x 0.25 mm ID x 0.25 µm df column.**

| Name           | CAS #       | R.T. (min) | R <sup>2</sup> | LOD  | Units |
|----------------|-------------|------------|----------------|------|-------|
| Etridiazole    | 2593-15-9   | 2.818      | 0.9994         | 0.50 | ng/g  |
| Atrazine       | 1912-24-9   | 3.750      | 0.9992         | 0.20 | ng/g  |
| Terbutylazine  | 5915-41-3   | 3.843      | 0.9996         | 0.20 | ng/g  |
| Terbacil       | 5902-51-2   | 3.963      | 0.9995         | 0.50 | ng/g  |
| Vinclozoline   | 50471-44-8  | 4.225      | 0.9992         | 0.20 | ng/g  |
| MGK 264        | 113-48-4    | 4.635      | 0.9992         | 0.20 | ng/g  |
| MGK 264 isomer | 113-48-4    | 4.708      | 0.9993         | 5.00 | ng/g  |
| Penconazole    | 66246-88-6  | 4.748      | 0.9999         | 0.20 | ng/g  |
| Fipronil       | 120068-37-3 | 4.793      | 0.9997         | 0.10 | ng/g  |
| Procymidone    | 32809-16-8  | 4.853      | 0.9998         | 0.20 | ng/g  |
| Pacllobutrazol | 76738-62-0  | 4.935      | 0.9994         | 0.20 | ng/g  |
| Flutriafol     | 76674-21-0  | 5.012      | 0.9996         | 0.50 | ng/g  |
| Fludioxonil    | 131341-86-1 | 5.100      | 0.9997         | 0.20 | ng/g  |
| Tricyclazole   | 41814-78-2  | 5.118      | 0.9991         | 1.00 | ng/g  |
| Myclobutanil   | 88671-89-0  | 5.165      | 0.9998         | 0.10 | ng/g  |
| Flusilazole    | 85509-19-9  | 5.183      | 0.9996         | 0.20 | ng/g  |
| Bupirimate     | 41483-43-6  | 5.198      | 0.9994         | 0.10 | ng/g  |
| Lenacil        | 2164-08-1   | 5.602      | 0.9997         | 0.50 | ng/g  |
| Hexazinone     | 51235-04-2  | 5.675      | 0.9998         | 0.20 | ng/g  |
| Tebuconazole   | 107534-96-3 | 5.693      | 0.9998         | 0.50 | ng/g  |
| Propargite     | 2312-35-8   | 5.713      | 0.9991         | 2.00 | ng/g  |
| Iprodione      | 36734-19-7  | 5.848      | 0.9994         | 0.20 | ng/g  |
| Pyriproxyfen   | 95737-68-1  | 6.146      | 0.9998         | 0.50 | ng/g  |
| Fenarimol      | 60168-88-9  | 6.312      | 0.9998         | 0.20 | ng/g  |
| Etofenprox     | 80844-07-1  | 6.885      | 0.9996         | 0.20 | ng/g  |
| Fluridone      | 59756-60-4  | 6.998      | 0.9992         | 0.10 | ng/g  |


Note: Captan, captafol, chlorfenapyr, folpet, cyprodinil, pyrimethanil, triadimefon, triadimenol, and triflumizole were incurred in the strawberry matrix or had chemical interferences that effected the low end of their calibration curves, so their linearity and LOD were not reported even though they were spiked and detected.

## Non-Target Full Mass Spectra All the Time

Analysis of 20 pg on column of bupirimate in strawberry on a 15 m x 0.25 mm ID x 0.25 µm df Rxi-5MS column demonstrating the software's ability to accurately identify analytes at low concentration in matrix.




## Examples of incurred pesticides in strawberry identified using NTD.



## NonTarget Deconvolution

Non-target analysis of a strawberry QuEChERS extract without cleanup and with instrument performance standards spiked at 500 pg. TIC and AIC highlighting representative analytes with a NIST library similarity score > 900.



Representative list of putatively identified peaks from the above figure.

| Name                         | CAS        | Formula   | R.T. (min) | Similarity | Type       | S/N  |
|------------------------------|------------|---|------------|------------|------------|------|
| N-Nitrosodimethylamine       | 62-75-9    | C <sub>2</sub> H <sub>6</sub> N <sub>2</sub> O                  | 1.949      | 958        | Non-Target | 86   |
| p-Methylcatechol             | 452-86-8   | C <sub>7</sub> H <sub>6</sub> O <sub>2</sub>                    | 2.301      | 949        | Non-Target | 1006 |
| Salicylic acid               | 69-72-7    | C <sub>7</sub> H <sub>6</sub> O <sub>3</sub>                    | 2.322      | 932        | Non-Target | 1138 |
| trans-Cinnamic acid          | 140-10-3   | C <sub>9</sub> H <sub>8</sub> O <sub>2</sub>                    | 2.705      | 937        | Non-Target | 1143 |
| 2,5-Di-tert-butylphenol      | 5875-45-6  | C <sub>14</sub> H <sub>22</sub> O                               | 2.962      | 918        | Non-Target | 4461 |
| Nerolidol                    | 7212-44-4  | C <sub>15</sub> H <sub>26</sub> O                               | 3.130      | 930        | Non-Target | 633  |
| Indole-3-acetaldehyde        | 2591-98-2  | C <sub>10</sub> H <sub>8</sub> NO                               | 3.611      | 901        | Non-Target | 1265 |
| PCB 18                       | 38444-73-4 | C <sub>12</sub> H <sub>7</sub> Cl <sub>5</sub>                  | 3.890      | 933        | Spiked     | 1518 |
| Pyrimethanil                 | 53112-28-0 | C <sub>12</sub> H <sub>13</sub> N <sub>3</sub>                  | 3.909      | 947        | Incurred   | 2668 |
| PCB 28                       | 55712-37-3 | C <sub>12</sub> H <sub>7</sub> Cl <sub>5</sub>                  | 4.187      | 937        | Spiked     | 1827 |
| n-Hexadecanoic acid          | 57-10-3    | C <sub>16</sub> H <sub>32</sub> O <sub>2</sub>                  | 4.384      | 946        | Non-Target | 808  |
| PCB 52                       | 41464-47-5 | C <sub>12</sub> H <sub>6</sub> Cl <sub>4</sub>                  | 4.404      | 919        | Spiked     | 1686 |
| Triphenylmethane             | 519-73-3   | C <sub>19</sub> H <sub>15</sub>                                 | 4.601      | 915        | Spiked     | 862  |
| Cyclic octatomic sulfur      | 10544-50-0 | S <sub>8</sub>  | 4.745      | 910        | Non-Target | 98   |
| Octadecanoic acid            | 57-11-4    | C <sub>18</sub> H <sub>36</sub> O <sub>2</sub>                  | 5.001      | 953        | Non-Target | 1069 |
| TDCPP                        | 13674-87-8 | C <sub>9</sub> H <sub>15</sub> Cl <sub>3</sub> O <sub>2</sub> P | 5.544      | 944        | Spiked     | 940  |
| 9-Octadecenamide, (Z)-       | 301-02-0   | C <sub>18</sub> H <sub>33</sub> NO                              | 5.589      | 905        | Non-Target | 1382 |
| Triphenyl phosphate          | 115-86-6   | C <sub>18</sub> H <sub>15</sub> O <sub>4</sub> P                | 5.729      | 915        | Spiked     | 837  |
| Stearic acid 1-monoglyceride | 123-94-4   | C <sub>21</sub> H <sub>42</sub> O <sub>4</sub>                  | 6.498      | 909        | Non-Target | 283  |
| Supraene                     | 7683-64-9  | C <sub>30</sub> H <sub>50</sub>                                 | 6.772      | 903        | Non-Target | 290  |
| Vitamin E                    | 59-02-9    | C <sub>29</sub> H <sub>50</sub> O <sub>2</sub>                  | 7.479      | 929        | Non-Target | 2721 |
| Campesterol                  | 474-62-4   | C <sub>28</sub> H <sub>48</sub> O                               | 7.718      | 914        | Non-Target | 549  |
| Sitosterol                   | 83-46-5    | C <sub>27</sub> H <sub>50</sub> O                               | 7.911      | 911        | Non-Target | 1611 |
| Isofucosterol                | 481-14-1   | C <sub>29</sub> H <sub>48</sub> O                               | 7.944      | 924        | Non-Target | 2329 |

Qualitative analysis may be used to screen for additional pesticides or to tentatively identify and characterize other aspects of the sample (e.g. in this case, flavor and fragrance characteristics of strawberries).

## Conclusions

- The Pegasus BT delivers a superior combination of quantitative and qualitative information in the same sample injection without sacrificing sensitivity.
- Target Analyte Find was demonstrated to quantitate across more than 4 orders of magnitude in food matrix, and up to 100 times lower than the regulatory guideline of 10 ng/g for the majority of pesticides in this study, with a run time of less than 9 minutes.
- LECO's industry-leading spectral deconvolution software provides uncontaminated mass spectra with unsurpassed spectral fidelity for library searching.

