### Introduction

Coffee is one of the most consumed beverages in the world and the industry around it is an important part of the global economy. As expected with commodities, there is a large amount of taste and flavor variation in coffee that can relate to differences in the variety and geographical origin of the beans, storage and processing conditions, roasting conditions, and brewing methods. An understanding of these differences can be helpful for quality control, process optimization, and also for providing information on flavors and characteristics that direct consumers to their preferred styles. In addition to the expected variation, the aroma profile for coffee is quite complex and comprised of a large number of individual analytes, creating an analytical challenge. Non-targeted chemical analysis techniques, like gas chromatography with mass spectrometry (GC-MS) and headspace solid phase micro-extraction (HS-SPME), are well-suited to address these challenges. With these methods, volatile and semi-volatile analytes were collected from the coffee samples, separated, and detected, resulting in identification and relative quantification information for hundreds of analytes. Analytes of interest do not need to be determined prior to acquisition, so the data were generally characterized to investigate the samples and their differences. Comprehensive two-dimensional gas chromatography ( $GC \times GC$ ) increases peak capacity and enhances S/N compared to GC, and also creates structured chromatograms. These additional analytical capabilities were explored and led to the detection of more analytes and an improved understanding of these complex samples. In this work, coffee brewed from a variety of beans was compared to investigate variations related to roast level.

### Methods

Coffee brewed from six different types of beans was compared with HS-SPME and GC or GC×GC coupled to TOFMS. The beans were from four geographical origins (Peru, Costa Rica, Kona, and Colombia) with a medium roast style from all four, and a dark roasted style from Costa Rica and Kona. The coffee was prepared by coarsely grinding each bean and brewing by French Press. Four ounces of water (100 °C) were added per 1 Tbs. whole beans and the coffee was pressed after 4 min of steeping. For HS-SPME analysis, 4 mL of coffee were transferred to a 20 mL vial which was incubated for 5 min at 60 °C, and then extracted with a DVB/CAR/PDMS fiber (Supelco) for 5 min at the same temperature. The samples were subsequently analyzed by GC-TOFMS and GC×GC-TOFMS with instrument conditions listed in Table 1. Single dimension GC data were acquired with the GC×GC column configuration by simply turning the modulator off, which allowed for rapidly switching between acquisition modes. Data for an alkane standard was also acquired for retention index calculations.

Table 1. GC-TOFMS and GC×GC-TOFMS	(Pegasus <sup>®</sup> BT 4D) Conditions
-----------------------------------	---

GC	Agilent 7890 with LECO Dual Stage Quad Jet
	Modulator and L-PALS Autosampler
Injection	SPME, 3 min desorption, split 5:1 in 250 °C inlet
Carrier Gas	He @ 1.4 ml/min, Corrected Constant Flow
Columns	Rxi-5Sil MS, 30 m x 0.25 mm i.d. x 0.25 $\mu$ m coating
	Rxi-17SilMS, 0.3 m x 0.25 mm x 0.25 $\mu$ m coating (Restek)
Oven Program	3 min at 40 °C, ramped 10 °C/min to 250 °C, hold 5 min
	Secondary oven maintained +10 °C relative to primary
Modulation	1.2 s, modulator temperature maintained +15 °C relative
	to secondary oven
Transfer Line	250 °C with uncoated guard column
MS	LECO Pegasus BT
Ion Source	250 °C
Temp	
Mass Range	33-510 m/z
Acquisition Rate	10 spectra/s (GC) and 100 spectra/s (GC×GC )

The use of GC×GC for these complex samples provided several key benefits with examples shown in Figures 1-3. GC×GC can: Increase peak capacity

The increased peak capacity provides better separations for complex samples as coelutions in the first dimension can often be separated in the second dimension. Figure 1 shows an example where a single peak marker was determined in the GC data that was revealed to be two analytes in the  $GC \times GC$ data. These analytes coelute in the first dimension, and the mass spectrum derived from the single dimension GC separation is the combination of the spectra of the two analytes and has a poor library similarity score. The analytes are chromatographically separated in the GC×GC data and have library matches with similarity scores into the 800s. What was one unknown with GC was determined to be two knowns with  $GC \times GC$ : 2,3-dimethyl pyridine with coffee and caramel odor notes and the furan, 5-methyl-2(5H) furanone.



# Aroma Profile of Coffee with GC, GC×GC, and TOFMS

Elizabeth M. Humston-Fulmer and Joseph E. Binkley | LECO Corporation, St. Joseph, MI USA

## **GC×GC** Benefits

• Enhance S/N

Create structured chromatograms

separated in the second dimension, as shown here.

GC×GC with thermal modulation is also expected to provide an enhancement in the S/N, demonstrated in Figure 2. This enhancement comes from thermal focusing at the modulator, which sharpens and focuses peaks just prior to detection. In the GC data, a single peak was identified as 3-phenyl furan. This furan is known to occur in coffee, and had a S/N above the threshold with both GC and  $GC \times GC$ . The  $GC \times GC$  data revealed that a second analyte was also present that was below the S/N threshold in the GC separation, but above the threshold after thermal focusing. The analyte, identified 5-hydroxymethyl furfural, has caramel and buttery odor properties and was only detected with the GC×GC separation.



Figure 2. GC×GC with thermal modulation leads to an enhanced S/N because of thermal focusing at the modulator.

The structured nature of the chromatograms is demonstrated in Figure 3. The complementary nature of the stationary phases leads to chromatograms where compound classes tend to elute in bands across the GC×GC separation space In this figure, peak markers are color coded by compound class for a collection of representative analytes to help visualize these bands. This aspect of GC×GC allows for rapid characterization of the samples and visual determination of the types of analytes present.



Figure 3. The complementary nature of the stationary phases leads to compound classes eluting in structured bands in the GC×GC separation space. This helps with visual characterization and highlights many of the types of analytes that were observed in these data.



# **GC×GC** Benefits

# **Representative Samples**

Beans from four geographical origins (Peru, Costa Rica, Colombia, and Kona) were analyzed, with a medium roast from each, and a dark roast from Costa Rica (French Roast), and Kona (Dark Kona). Representative TIC chromatograms for each coffee sample are shown in Figure 4.



Figure 4. Representative TIC chromatograms of each coffee sample.

### Sample Differentiation

The coffee samples have many similarities, but some differences can also be noted. The structured nature of the plot gives insight to these differences even before peak finding has been performed. For example, the medium roast Peru sample has some unique peaks in the alkane band that are not present in the other samples. The dark roast Dark Kona sample has more intensity in the nitrogen-containing ring band compared to the other samples. Automated data processing provided peak information for the entire sample, and a wide range of compound types were found including alkanes, terpenes, aldehydes, ketones, furans, nitrogencontaining rings, aromatic compounds, and thiophenes. All identifications were tentatively determined with retention index and spectral similarity to library databases. Having information on hundreds of individual analytes allowed for making comparisons of the coffee prepared from the different beans. Peak areas were compared for each analyte across the sample set to observe general trends and differences. A collection of information for some of the sampledistinguishing analytes observed here is shown in Table 2. These differences appear to relate to roast level. Of note, several specific alkanes were observed elevated in the Peru samples and several specific pyridines were observed elevated in the Dark Kona samples, consistent with the preliminary observations based on the chromatographic structure. The aroma properties were also compiled for these sample-distinguishing analytes. Several analytes that were observed at higher levels in the coffee from darker roasted beans had odor descriptors like roasted, coffee, smoke, and burnt; and several analytes at higher levels in the coffee from the medium roasted beans had odor descriptors like caramel, buttery, baked bread, and nutty. Many other analytes were observed that would also contribute to the overall aroma profile.

Table 2. Some analyte differences that appear to relate to roast level.

	tR	Similarity	Formula	CAS	RI (Obs)	RI (LIb)	Peru	Costa Rica	Kona	Colombian	French Roast	Dark Kona	Odor and flavor notes
decane, 3,7-dimethyl-	673.2	867	C12H26	17312-54-8	1124	1125							
undecane, 2-methyl-	710.4	908	C12H26	7045-71-8	1163.8	1164							
undecane, 3-methyl-	716.4	949	C12H26	1002-43-3	1170	1170							
2,3-pentanedione	195.0	929	C <sub>5</sub> H <sub>8</sub> O <sub>2</sub>	600-14-6	707.3	698							caramel, nutty, sweet, butter, creamy, cheese, pungent
furfural	356.0	966	C <sub>5</sub> H <sub>4</sub> O <sub>2</sub>	98-01-1	835.9	833	•						woody, almond, fragrant, baked bread, sweet
2(5H)-furanone	454.6	896	$C_4H_4O_2$	497-23-4	916.6	918							buttery
2-furancarboxaldehyde, 5-methyl-	510.8	950	C <sub>6</sub> H <sub>6</sub> O <sub>2</sub>	620-02-0	966.2	965							spice, caramel, maple
3(2H)-furanone, dihydro-2-methyl-	323.5	938	C <sub>5</sub> H <sub>8</sub> O <sub>2</sub>	3188-00-9	809.8	809					_		bread, buttery, nutty, sweet, solvent
thiazole, 2-methyl-	323.5	839	C <sub>4</sub> H <sub>5</sub> NS	3581-87-1	809.8	815							green, vegetable
furan, 2,2'-methylenebis-	637.9	935	C <sub>9</sub> H <sub>8</sub> O <sub>2</sub>	1197-40-6	1086.7	1088				_			rich, roasted
acetophenone	621.2	894	C <sub>8</sub> H <sub>8</sub> O	98-86-2	1070.3	1065							sweet, almond, pungent, hawthorn, mimosa, acacia, chemical
pentanoicacid, 4-oxo-, methyl ester	537.2	912	$C_6H_{10}O_3$	624-45-3	989.4	982							caramellic
furan, 2-(2-furanylmethyl)-5-methyl-	728.4	926	C10H10O2	13678-51-8	1183	1190							
2(3H)-furanone, dihydro-5-methyl-	500.1	899	C <sub>5</sub> H <sub>8</sub> O,	108-29-2	956.7	958							cocoa, woody, sweet, herbal, warm, tobacco
1-(2-thienyl)-1-propanone	736.4	854	C <sub>7</sub> H <sub>8</sub> OS	13679-75-9	1190.7	1185				_	-		caramel, creamy
furan, 2,2'-[oxybis(methylene)]bis-	836.0	915	C10H10O3	4437-22-3	1305.5	1299							coffee, nutty, earthy
benzoxazole, 2-methyl-	670.4	839	C <sub>8</sub> H <sub>7</sub> NO	95-21-6	1120.1								burnt, tobacco, phenolic, meaty, powdery, capers
1H-pyrrole, 1-ethyl-	331.2	923	C <sub>6</sub> H <sub>9</sub> N	617-92-5	816.6	821							burnt flavors
methyl 2-furoate	524.0	888	C <sub>6</sub> H <sub>6</sub> O <sub>3</sub>	611-13-2	977.8	980							tobacco, fruity, mushroom, fungus, sweet
phenol, 2-methoxy-	644.0	941	C <sub>7</sub> H <sub>8</sub> O <sub>2</sub>	90-05-1	1092.6	1090							smoke, spice, vanilla, woody, phenolic
5,6,7,8-tetrahydroquinoxaline	759.2	845	$C_8H_{10}N_2$	34413-35-9	1216.4	1223							roasted, nut, musty, bean, cereal, corn, chip, cheese
phenol	529.9	949	C <sub>6</sub> H <sub>6</sub> O	108-95-2	983.1	980							phenolic plastic rubber
methyl-2-thiophene carboxylate	665.6	852	$C_6H_6O_2S$	5380-42-7	1115								burnt flavors
dihydro-2(3H)-thiophenone	554.1	922	C₄H <sub>6</sub> OS	1003-10-7	1004.7								burnt, garlic
phenol, 4-ethyl-2-methoxy-	818.0	900	$C_9H_{12}O_2$	2785-89-9	1284.3	1282							spicy, smoky, bacon, phenolic, clove
pyridine, 2-ethyl-	444.6	853	C <sub>7</sub> H <sub>9</sub> N	100-71-0	908.2	906							green, grassy
pyridine, 3-methoxy-	551.6	831	C <sub>6</sub> H <sub>7</sub> NO	7295-76-3	1002.3	1005							
pyridine, 3-methyl-	391.9	843	C <sub>6</sub> H <sub>7</sub> N	108-99-6	864.9	863							earthy, hazelnut, green flavors
furfuryl sulfide	969.3	888	C <sub>10</sub> H <sub>10</sub> O <sub>2</sub> S	13678-67-6	1474.6	1463							coffee, mushroom, earthy, powerful, meaty, sulfury
pyridine, 3-ethyl-	503.5	912	C <sub>7</sub> H <sub>9</sub> N	536-78-7	959.9	959	0.5	1 1.5	2 2.5	3 3	3.5 4 4	1.5 5	roasted flavors; tobacco, oakmoss, leather odors

### Conclusions

The Pegasus BT 4D is a powerful analytical tool that allows for non-target GC or GC×GC-TOFMS analyses to gain insight and learn more about your complex samples. In this work, variations in coffee samples relating to the roast level of the beans were investigated. HS-SPME collected the volatile and semi-volatile analytes from the coffee samples, GC or GC×GC separated the analytes from each other, and TOFMS detection provided identification and relative quantitation information for hundreds of analytes. A comparison of six different coffee samples prepared from medium and dark roasted beans found specific analytes that appear to relate to the roasting styles, independent of geographical origin. The benefits of GC×GC were demonstrated and provided information about these complex samples, and uncovered some analytes that were difficult to measure with GC alone.