



Delivering the Right Results

# Methods for Profiling Gin for Essential Oil Components with GC-TOFMS and GCxGC-TOFMS

Donald C. Hilton, Megan McGuigan, and Scott Pugh • LECO Corporation, St. Joseph, Michigan

## INTRODUCTION

The detection of minor components in gin is important in the evaluation of the flavor composition and/or the geographic origin of the product. In this investigation we demonstrate the use of a single injection for the identification of flavor components in a set of six gin samples. Neat samples were examined by injection without further preparation and analyzed by both gas chromatography time-of-flight mass spectrometry (GC-TOFMS) and comprehensive two-dimensional gas chromatography time-of-flight mass spectrometry (GCxGC-TOFMS). Analysis by GC-TOFMS has the advantage of being simple, fast, and sensitive to a full range of analytes. Analysis by GCxGC-TOFMS allows for the detection of trace level analytes that may not be seen using other methods. Chromatographic peaks generated in high-speed GC and GCxGC separations are often quite narrow and require the use of a relatively fast detector in order to be fully characterized. The high acquisition speed of the TOFMS allows for the characterization of these narrow peaks and also allows for spectral deconvolution of overlapping chromatographic peaks. The mass spectra were automatically compared against the NIST and Terpene<sup>1</sup> libraries for peak identification.

## SAMPLE

By definition, all gin samples are derived from the juniper berry<sup>2</sup>. Gin is made through the distillation of the juniper berry or by adding juniper berry oil to a refined spirit. Many gins also contain additional essential oils such as citrus oils in order to enhance their flavor. In this study, a set of six gin samples was studied and compared to a juniper oil standard. Initial data review was done on a GC-TOFMS analysis of juniper oil. The identified compounds in the juniper oil were then used to build a reference against which all other gin samples were compared. This was done through the use of the "Reference" feature in ChromaTOF<sup>®</sup> which allows for the automated comparison of chromatograms against a standard based on retention time, a user-defined s/n threshold and spectral similarity. Once the comparison was made on the one-dimensional data, the two-dimensional chromatograms were compared to their one-dimensional counterparts. The one-dimensional data shows a good correlation between the juniper oil standard and the gin samples. The two-dimensional chromatograms show enhanced selectivity and separation of target compounds in addition to increased s/n for target analytes and an overall increase in the total number of compounds identified.

## EXPERIMENTAL

### GC-TOFMS Analysis

#### Sample Introduction:

Agilent Autosampler with a 5 L syringe  
Injection: 1 µL (splitless)

#### GC: Agilent 6890 Gas Chromatograph

Inlet: split/splitless; 275°C  
Carrier Gas: He at 1.5 mL/min  
Column: 10 m x 0.18 mm x 0.20 m RTX-5  
GC Oven: 60°C (1 min hold) programmed to 250°C at 20°C/min (2 min hold)  
MS Transfer Line: 250°C

#### MS: LECO TruTOF<sup>®</sup> HT

Ionization: EI at -70 eV  
Ion Source Temperature: 300°C  
Spectral Acquisition Rate: 10 spectra/s  
Acquired Mass Range: 35-500  
Acquisition Delay: 60 s

#### Instrument Control and Data Review:

ChromaTOF optimized for TruTOF HT

### GCxGC-TOFMS Analysis

#### Sample Introduction:

Agilent Autosampler with a 10 L syringe  
Injection: 1 µL (splitless)

#### GC: Agilent 6890 Gas Chromatograph

Inlet: split/splitless; 275°C  
Carrier Gas: He at 0.6 mL/min  
Primary Column: 10 m x 0.18 mm x 0.18 m DB-5  
Secondary Column: 1 m x 0.10 mm x 0.10 m BPX-50  
Primary Oven: 80°C (2 min hold) programmed to 225°C at 5°C/min (3 min hold)  
Secondary Oven: 5°C offset from the primary oven  
Thermal Modulator: 30°C offset from the primary oven  
Modulation Period: 5 s  
MS Transfer Line: 250°C

#### MS: LECO Pegasus<sup>®</sup> 4D

Ionization: EI at -70 eV  
Ion Source Temperature: 200°C  
Spectral Acquisition Rate: 150 spectra/s  
Acquired Mass Range: 35-350  
Acquisition Delay: 100 s

#### Instrument Control and Data Review:

ChromaTOF optimized for Pegasus 4D



## RESULTS AND DISCUSSION




Figure 1. Total Ion Chromatogram (TIC) and peak table for the analysis of Juniper Oil.




Figure 2. Total Ion Chromatogram (TIC) for the analysis of the Gin\_A sample by GC-TOFMS. The inset chromatogram shows the detail for the smaller peaks along with the (a) Peak True and (b) Reference Spectra for one of the compounds, Caryophyllene.

Figure 2 shows a TIC for the analysis of a commercially-available gin sample. The large peak at the beginning of the chromatogram is Ethanol. The inset shows the smaller peaks in detail. Shown above the chromatogram are the (a) Peak True and (b) Reference Spectra for one of the analytes, Caryophyllene. The Peak True spectrum contains all of the ions associated with a particular compound following the application of the Deconvolution algorithm. The peak true spectrum was then compared against the reference standard in order to identify the analyte as Caryophyllene. In this case, the Reference Spectrum is the Caryophyllene Peak True spectrum from the Juniper Oil chromatogram shown in Figure 1.

A total of 17 compounds that were present in the Juniper Oil standard were identified in the "Gin\_A" standard. An additional 43 compounds were identified with a s/n greater than 25. Of those compounds, 22 were identified following comparison with the standard NIST and Terpene libraries.

### Observations

Figure 1 shows the TIC for a GC-TOFMS analysis of a juniper oil standard that was completed in just under 500 s. Note that this sample was diluted 1:100 in MeOH in order to prevent sample overload. Automated data processing revealed 115 compounds with a s/n of 25 or greater. All spectra were automatically compared against the NIST and Terpene libraries for peak identification. Shown in the table to the left is the peak table for all analytes with a library match of 750 or greater. A total of 76 compounds were identified and these were used to build a reference, against which all of the 1-D gin data was compared. The reference was set with constraints allowing for a ±0.5 s shift in retention time, a s/n threshold of 10 and a spectral similarity of at least 600. The data processing method for each of the gin samples analyzed by GC-TOFMS included this reference for peak identification. Additional compounds not present in the reference standard were identified via comparison with the NIST and Terpene libraries.

## SPECTRAL DECONVOLUTION




Figure 4. Example of spectral deconvolution for the identification of two overlapping components in the Juniper Oil sample.

### Observations

Figure 4 shows an example of spectral deconvolution observed in the analysis of the juniper oil standard. The TIC (red) appears to show the presence of a single chromatographic peak. However, by plotting the extracted ions (m/z 93 and 177) it is clear that there are in fact two analytes present. The spectra on either side of the chromatogram are for the two overlapping analytes. The Caliper spectra show all ions present at the peak apex. As a result of the overlap between the two compounds, ions for both are present in each caliper spectrum. Following deconvolution, the Peak True spectrum is generated which contains ions specific to that compound. Note the presence of m/z 121 and 136 in the caliper spectrum on the right. These are both due to the overlap between the two compounds. Following deconvolution, the Peak True spectrum on the right no longer contains these ions. The Peak True spectra were automatically compared against library spectra in order to identify the compounds as 2-Carene and (2-methyl-2-propenyl)-Benzene with similarity rankings of 910 and 775 respectively.

## COMPARISON OF SAMPLES

Compound	m/z	Gin_A GC	Gin_C GCxGC	Gin_B GCxGC	Gin_C GCxGC	Gin_D GCxGC	Gin_E GCxGC	Gin_F GCxGC						
α-Cubebene	161	10.37	399.44	17.14	36.00									
Cyclohexane, 1-ethenyl-1-methyl-2,4-bis(1-methylheptyl)-	81	10.37	399.44	11.82	70.47	15.00	319.34	10.18	228.51					
Isolindene	161	28.45	144.54			30.63	59.94	20.68	38.06					
Caryophyllene	133	91.00	716.59	18.12	98.69	137.75	272.97	179.65	427.18					
β-Element	121	109.14	808.41	36.92	216.14	249.96	742.49							
α-Caryophyllene	93	150.65	2166.20	42.59	510.62	212.81	1311.90	295.39	1416.00					
(E)-Epi-bicyclosequiphellandrene	161							40.33	51.51					
Germacrene D	161	148.36	669.50	41.52	165.02	424.57	865.34	360.99	718.37					
β-Selinene	93					12.03	86.17							
α-Muurolene	105	44.36	572.71			60.52	305.16	56.51	247.48					
d-Cadinene	161	78.35	281.76			88.45	152.00							
β-Cadinene	161	368.27	1170.50	66.24	177.94	482.03	626.95	135.52	256.27	305.23	328.56			
α-Cadinene	119	11.28	82.44			12.97	168.23							
Caryophyllene oxide	79	16.84	96.99											
α-Cadinol	93	39.67	107.78			31.39	47.73	86.23	238.26					
β-Bisabolene	69	21.56	110.78	20.85	77.60			21.38	104.80					
Limonene	68	539.61	2344.8	104.84		11.08		208.7	3131.3	38.59	809.2			
Linalool	71	1142.6	13627	540.97	3930	192.8	1651.1	2057.5	12297	2413.5	17469	1843.9	13160	
Verbenyl Ethyl Ether	100	76.801	545.88	53.287	234.17	3		36.711	61.912	43.3	175.12	21.275		
Limonene	95	56.698	642.31	29.238	192.34			92.006	141.91	632.86	149.33	732.59	148.26	718.84
Total # of peaks detected		60	462	34	356	27	300	63	325	129	771	67	421	

Table 1. Comparison of the s/n values observed for a select group of compounds in 6 different gin samples analyzed by GC-TOFMS and GCxGC-TOFMS.

### Observations

Table 1 shows a comparison of each of the six gin samples analyzed by GC-TOFMS and GCxGC-TOFMS. The top part of the table shows a select group of compounds that were identified based on the reference comparison to juniper oil. The bottom set of compounds were not in the juniper oil, but seemed to be common among the gin samples and are likely flavor additives. As expected, there was a significant increase in the s/n values observed using GCxGC-TOFMS. Gin sample C did not have any similarities to the juniper oil standard. Gin sample E did not have many similarities to the juniper oil standard either, but contained many other essential oil components which was expected as this was labeled as a botanical gin.

## CONCLUSIONS

This study showed a comparison of GC and GCxGC experiments for the screening of gin samples for their flavor components. Liquid injections were done of the neat samples with no sample preparation required. Automated data processing revealed a large number of compounds in samples analyzed by both methods. The use of the Reference feature in ChromaTOF allows for the fast and easy comparison of a group of samples to a reference sample, which in this case was the Juniper Oil. The fast detection offered by the TOFMS makes spectral deconvolution possible and allows for the detection and identification of analytes that may not be chromatographically resolved. As expected, a significant increase in s/n was observed in the GCxGC results relative to the GC results; however, the GC analysis revealed a large amount of information about the samples as well.

<sup>1</sup> The Terpene Library contains mass spectra of essential oil components compiled by Robert P. Adams, Baylor University Plant Biotechnology Center.  
<sup>2</sup> K.H. Newman et al. Chromatogr. A 1164 (2002) 351-370