

Analysis of Synthetic Drugs Using Electron and Chemical Ionization High Resolution Time-of-Flight Mass Spectrometry

Introduction

The development of unregulated recreational drugs has been increasing at an alarming rate since the mid-2000s. Relatively simple organic transformations produce novel and licit psychotics that can elude detection by standard analytical methods (e.g., GC-MS). This is a major problem for law enforcement agencies, forensic laboratories, and the medical community since the legal status, physiological effects, and long-term toxicity of these synthetic drugs are unclear. Detection and characterization of these substances through routine forensic analyses is challenging due to their novelty and lack of appropriate instrumental methods. Matters are complicated further by the wide range of active ingredients and variety of botanical blends in samples which can be difficult to handle instrumentally.

High performance time-of-flight mass spectrometry (TOFMS) is a practical choice for the analysis of these substances:

- Data acquisition is comprehensive and robust
- Single acquisition, high-resolution accurate mass data can be probed multiple times for traditional and emerging drugs
- High resolving power minimizes background interferences
- Excellent mass accuracy values facilitate confident elemental formula determination for compound characterization
- Complete and high-quality spectral data can be searched against databases

Traditional



Synthetic



Figure 1: Traditional and Synthetic Drug Examples

In this study, two samples (Cases 1 and 2) obtained from a collaborating forensics laboratory were analyzed and found to contain traditional and new psychotropic substances. The combination of electron impact and chemical ionization high resolution TOFMS was vital for unambiguous identification of new drugs in these mixtures.



GC Parameters

GC: Agilent 7890 and 7693 Auto Sampler Column Type: Restek Rxi-5Sil MS (30 m, 0.25 mm ID, 0.25 mm df) Injection: 0.5 mL, Split 150:1 (CI 100:1); Inlet Temp. 250°C Oven: 60°C (2 min) to 300°C at 20°C/min (1 min) **Carrier Gas:** He, Constant Flow (1.00 mL/min)

MS Parameters

Spectrometer: Ion Sources: Folded Flight Path: Spectral Acquisition Mass Range (m/z): Calibration: CI Reagent Gas:

Analysis Workflow





David E. Alonso, Joe Binkley, and Jeff Patrick | LECO Corporation, Saint Joseph, MI USA

Experimental

Sample Preparation

Case 2: Drug Pipe (Cannabinoids) Dark Mixture (93 mg) 3 mL 2:1 DCM/MeOH N-C₅H₁₁ N-C₅H₁₁

Figure 2. Sample Preparation Procedures

Instrument Parameters

LECO Pegasus[®] GC-HRT LECO EI, CI HR (R = 25,000 FWHM) 6 spectra/second 30-650 (CI 60-650) **PFTBA** (Internal) 5% Ammonia in Methane

Figure 3: Drug Analysis Workflow; EI and CI-HRT Analysis

Analysis of the sample produced the analytical ion chromatogram shown in Figure 4. The four major components of the mixture identified through library searches were 2-Fluorophenyl 2-methylbenzoate, N-Ethylamphetamine, Ethylone and Dimethylone (Figure 4). Library matches for these compounds were very good (799 to 935/1000); however, quick inspection of the CI-HRT data demonstrated that the number one hit for peak #2, N-ethylamphetamine ($C_{11}H_{17}N$), was incorrect (Figure 5). A formula search on the analyte's protonated molecular ion at m/z = 192.13825 resulted in the formula $C_{12}H_{18}NO$ (-0.21 ppm). An internet search using the formula $C_{12}H_{17}NO$, led to Cayman Chemical and Scientific Working Group Sites and 4-methylethcathinone as a potential candidate for this compound.





El and CI-HRT data for ethylone and dimethylone is shown in Figs. 6 and 7. Protonated molecular ions $(C_{12}H_{15}NO_3)$ for these structural isomers were found at m/z = 222.11249 (0.08 ppm) and m/z = 222.11247 (-0.015 ppm) respectively.



Figure 6: El (top) and CI-HRT Data (bottom) for Ethylone, Sample 1

Results (Case 1: White Powder)



Figure 5: EI (top) and CI-HRT Data (bottom) for Analyte 2, Sample 1

Figure 7: EI (top) and CI-HRT Data (bottom) for Dimethylone, Sample 1

Sample 2 was a heterogeneous, dark, solid material that was scraped from the inside of a confiscated drug pipe. A representative set of the over 130 compounds in the sample are listed in Table 1. These include nitriles, phenols, sterols and various heterocyclic compounds. An extracted ion chromatogram XIC (Figure 8) of the mixture shows the major component A, which was not present in any of the commercial libraries used during data processing. The XIC also displays an additional unknown (*B), as well as various drugs detected in the sample.



CI-HRT data was critical for the identification of unknowns A and B (Figure 9). EI-HRT Library hits for these unknowns were poor. Confirmation of unidentified compound formulas was accomplished using CI-HRT data. Internet and database searches of the data sets led to AM-2201 and JWH-022 as probable hits, which were then confirmed through the analysis of standard spectra.

A combination of EI and CI-HRT resulted in confident identification of traditional and emerging synthetic drugs. The GC-HRT's high resolution, quality spectral data, and excellent mass accuracy values are very useful for the identification of novel compounds in complex matrices. The accurate mass data provided formulas that were searched against forensic databases for identification of unknowns in drug samples.



0.00066

Results (Case 2: Drug Pipe)

Table 1: Compounds In Sample 2

Ave. = 0.53 ppm



Figure 9: EI (top) and CI-HRT Data (bottom) for Unknowns A and B

Summary