# Characterizing Extractables from Common Pharmaceutical Packaging Material by High Resolution Time-of-Flight Mass Spectrometry and Enhanced Gas Chromatography Separations

## INTRODUCTION

The characterization of extractable and leachable components from a wide GCxGC improves the peak capacity of a separation, which allows for Incorporating EI and CI accurate m/z information from HR-MS with GCxGC, variety of samples were evaluated with GCxGC and HR-TOFMS. range of materials is an important area of research. Information about chromatographically isolating more analytes from each other and from which generally provides cleaner spectra for interpretation, can yield better Representative samples and analytes are shown. extractable and leachable components from packaging and delivery devices interferences. This can help uncover new analyte peaks, provide spectra with identifications. for pharmaceutical products is a particular area of growing interest. Analytical fewer interferences, and lead to more identified analytes. Figure 6. GCxGC-HR-TOFMS data for five different testing results on this topic are part of many regulatory submission requirements samples are shown. One plastic syringe without a rubber stopper (A), two syringes with rubber stoppers (B and C to the FDA. USP 1663 provides guidance on extractables testing and a variety CI spectrun Plastic Syringe with Rubber Stopper and two rubber stoppers (D and E) were evaluated Observed Calc Mass Mass Formula MA ppm of analytical approaches can meet compliance. High resolution MS is often GCxGC was important for separating the com these complex samples, and HR-TOFMS (both El and CI) considered necessary for identification of unknowns, and sample complexity was necessary for the identifications. Several analyte and low-level detection continue to challenge these analyses. Here, we examples are shown in the associated table. demonstrate a workflow that uses comprehensive two-dimensional gas chromatography (GCxGC) with HR-TOFMS to help address these challenges.

### METHOD

Representative extract samples were prepared from various materials that are commonly used in packages and closures for pharmaceutical products. Butyl rubber stoppers and plastic syringes (with and without rubber components) were extracted with methylene chloride at room temperature for 72 h. The extracts were analyzed by GCxGC-HR-TOFMS (Pegasus<sup>®</sup> HRT<sup>+</sup> 4D, LECO), as described in Table 1. An alkane standard was also analyzed for retention index (RI) determinations. EI and CI (methane) data were collected with a multimode ionization source and used to support identifications.

Table 1. Instrument Conditions	
AS	LECO L-PAL3 Autosampler
Injection	1μL
GCxGC	LECO GCxGC QuadJet <sup>TM</sup> Thermal Modulator
Inlet	280 °C, splitless
Carrier Gas	He @ 1.40 mL/min, constant flow
Columns	Column 1: Rxi-5ms, 30 m x 0.25 mm i.d. x 0.25 µm coating (Restek) Column 2: Rxi-17Sil MS, 0.9 m x 0.25 mm i.d. x 0.25 µm coating (Restek)
Temperature Program	2 min 50 °C, ramp 8 °C/min to 340 °C, hold 5 min Secondary Oven: + 20 °C
Modulation	3 s with temperature maintained +15 °C relative to 2nd oven
Transfer Line	350 °C
MS	LECO Pegasus HRT <sup>+</sup> 4D
Source	LECO MMS <sup>TM</sup> (El and Cl, with methane)
Ion Source Temp	250 °C (El) and 165 °C (Cl)
Mass range	35-900 m/z (EI) and 60-900 (CI)
Acquisition Rate	125 spectra/s



Elizabeth M. Humston-Fulmer, David E. Alonso, and Joseph E. Binkley | LECO Corporation, Saint Joseph, MI USA

# ENHANCED SEPARATION WITH GCxGC




Figure 1. GC and GCxGC chromatograms for an extract of a plastic syringe indicate improved chromatographi separation. More analytes are isolated from each other in chromatographic space and separated from interferences. Additionally, structured chromatograms and elution order in the second dimension add support to identifications along with similarity score, RI in the first dimension, and accurate mass molecular formula and fragment determinations. As an example, Chemanox 22, a rubber antioxidant, can be observed separated from background interferences.





Figure 2. In some instances, GCxGC helps uncover new analytes that were hidden in a GC separation. Here, a phenol compound and dodecyl acrylate completely coelute and are combined as one peak in the 1D separation. With GCxGC, the analytes are chromatographically separated in the second dimension, and both were determined. This revealed more information than could be determined with just GC. Additionally, CI data added molecular ion support for dodecyl acrylate that did not have a molecular ion in the El data and also confirmed molecular ion for 2,4-di-t-butyl-6-nitrophenol.

# **IMPROVED IDENTIFICATIONS**



as shown here. HR-MS, both EI and CI, add crucial identification information. In instances where a molecular ion is no present in the EI data, the CI data can often support the formula determination, as shown here for a phthala observed in a plastic syringe




Figure 4. Accurate mass and formulae determinations can also allow for selecting an improved spectral match from the library database. In this case from a plastic syringe, both library matches are supported by nominal mass fragmentation patterns, but the formula for the first library hit was not well-supported and the second library hit was.




Figure 5. Accurate mass information can also help determine formulae for features that are not present in library databases. In this example, butyl oligomers that were not in NIST library databases were observed in a rubber stopper extract. Accurate mass information was used to determine the formula and elution in the structured space supported the identifications.





Figure 6. Chlorinated butyl oligomers that were not in NIST library databases were also observed in a rubber stopper extract. Accurate mass information was used to determine the formula and elution in the structured space supported the identifications.

In this work, GCxGC-HR-TOFMS was used to evaluate extracts from several pharmaceutically relevant materials. Two dimensions of separation helped address sample complexity, and high-resolution MS helped address analyte identification requirements. Several examples to highlight these benefits are shown, as well as representative samples and representative analytes.



#### SAMPLE CHARACTERIZATION



# CONCLUSIONS