

A Workflow for Comprehensive Analysis of Forensic Samples Using Electron and Chemical Ionization High Resolution Time-of-Flight Mass Spectrometry

David E. Alonso¹, Joe Binkley¹, and John Rorabeck² | ¹LECO Corporation, St. Joseph, MI; ²Berrien County Forensic Laboratory, Berrien Springs, MI

Introduction

Background

Controlled substances such as marijuana, cocaine, heroin and methamphetamines are commonly seized by law enforcement officials and routinely analyzed by forensic laboratories. More recently, synthetic drugs such as bath salts and cannabinoids, as well as naturally occurring "herbal remedies" and hallucinogenic mushrooms, have gained popularity as chemical stimulants and opiate substitutes. This study focuses on the comprehensive profiling of complex botanicals obtained from drug seizures.



Botanicals Bath Salts Magic Mushrooms

Analytical Challenges



- Chemical Diversity of Compounds
- Novelty of Synthetic Compounds
- Concentration Range of Sample Constituents
- Complexity of Botanical Matrices
- Inappropriate Sample Preparation Methods
- Unsuitable Instrumental Analysis Protocol
- Absence of Library Spectra for Unknowns
- Lack of Standards

The Solution

- **Tailored Sample Preparation**
- **GC-HR TOFMS – Pegasus® GC-HRT:**
 - High Quality Spectral Data
 - Comprehensive
 - Search Against Well-Established Databases (e.g., NIST, Wiley)
 - Excellent Mass Accuracy Values (<1 ppm) = Robust Formulas for Fragment, Molecular and Adduct Ions
 - High Resolution Deconvolution™ (HRD™)
 - High Resolving Power (up to 50,000)



EI/CI-HRT Workflow:



Experimental

Sample Preparation Methods

Botanical samples were obtained from a forensic laboratory after they had been analyzed and pertinent cases were completed. General analysis for volatile materials was accomplished through head-space solid phase microextraction (HS-SPME). Solvent extraction followed by derivatization was utilized for comprehensive profiling of samples.



HS-SPME (Volatiles)

Placed 0.05 g of crushed sample (mortar & pestle) in 10 mL headspace vial. Added 0.5 g NaCl and 1 mL of either deionized water or 10% NaOH (aq). Extraction was carried out using a 50/30 mm DVY/CAR/PDMS Stableflex, 24 Ga SPME fiber. The procedure included fiber conditioning (1 hr at 270°C), incubation (3 min at 95°C), extraction (15 min at 95°C, 250 rpm interval agitation (5s on, 2s off)), followed by desorption in the GC-injector (5 min at 250°C).

Solvent Extraction/Derivatization (Profiling)

Samples (0.025 g) were mixed with methanol (2 mL) and sonicated for 25 minutes. The heterogeneous mixture was filtered and dried with N₂ gas. The samples were then treated with MEOX (25 µL, 60°C, 1 hr), MSTFA (50 µL, 60°C, 1 hr) and transferred to 2 mL GC vials for analysis.


Instrument Parameters

GC	
Agilent 7890 with Gerstel MPS Auto Sampler	
Column	Restek Rxi-5 MS (30 m x 0.25 mm x 0.25 µm)
Carrier Gas, Flow	He, 1.0 mL/min Constant Flow
Injection	1 µL, Split 5:1, SPME
Inlet Temperature	270°C
Temp. Program	70°C (2 min) to 280°C at 20°C/min (3 min)


MS	
LECO Pegasus GC-HRT	
Transfer Line Temp.	300°C
Ion Source Temp.	EI 250°C; CI 200°C
Ionization	EI (70 eV); CI (140 eV)
Mass Range	EI 35 to 510; CI 60 to 510
Acquisition Rate	6 spectra/second
Calibration (Internal)	PFTBA
Reagent Gas	5% Ammonia in Methane

HS-SPME Results



A) Marijuana (High Quality Spectral Data)



B) Compressed Hash (Excellent Mass Accuracy)




C) Khat (High Resolution Deconvolution)





Liquid Extraction/Derivatization Results

D) Pokeweed (Discovery & Confirmation)



E) Magic Mushrooms (Comprehensive Profiling)



Summary

- EI/CI workflow facilitates confident compound identification
- Proper sample preparation is critical for the analysis of complex samples
- GC-HRT provides high quality, accurate mass data for:
 - Spectral similarity searches (NIST, Wiley, etc.)
 - Formula determination (fragment, molecular, and adduct ions)

LECO and Pegasus are registered trademarks of LECO Corporation. High Resolution Deconvolution and HRD are trademarks of LECO Corporation.