

# Comparison of Smoker's and Non-Smoker's Urine Using Comprehensive Two-Dimensional Gas Chromatography High Performance Time-of-Flight Mass Spectrometry

David E. Alonso\* and Joseph E. Binkley | LECO Corporation, St. Joseph, Michigan U.S.A.

## Introduction

Urine is a favored biofluid for diagnostic testing (Urinalysis) because it is non-invasive and large volumes are easily obtained. In addition, urine is relatively free from interfering proteins and lipids, and it tends to "hold" high concentrations of drugs and metabolites over extended periods of time. Modern, routine clinical tests include the determination of specific gravity, measurement of glucose, nitrates, etc.

In this study, a novel analytical approach was utilized for the effective characterization of compounds in two standard reference materials, NIST smoker's and non-smoker's urine.

## Objectives

- Implement the use of enhanced, comprehensive two-dimensional gas chromatography (GCxGC) for the separation of compounds in urine
- Use a benchtop, high performance time-of-flight mass spectrometer and powerful processing software to quickly and confidently identify compounds in urine
- Use software tools to compare smoker's and non-smoker's urine

## Sample Preparation

- 600 µL of urine treated with urease (37 °C, 15 min)
- The mixture was vortexed (2 min) and then centrifuged (12,000 g for 10 min)
- 200 µL of supernatant was transferred to a 2mL GC vial and evaporated to dryness (Speed Vac)
- The dry material was derivatized using a two-step procedure
  - Methoximation (20 µL of MEOX, 80 °C, 30 min)
  - Silylation (75 µL MSTFA, 80 °C, 30 min)

## Technology




Table 1. Instrument acquisition parameters

| Gas Chromatograph                   | Agilent 7890, LECO Dual Stage Quad Jet Modulator & L-PAL 3 Autosampler  |
|-------------------------------------|---|
| Injection                           | 1µL, Split 20:1, 280°C  |
| Carrier Gas                         | He @ 1.4 ml/min, Constant Flow  |
| Columns (1 <sup>st</sup> Dimension) | Rxi-5 MS, 30 m x 0.25 mm i.d. x 0.25 µm (Restek, Bellefonte, PA, USA)   |
| (2 <sup>nd</sup> Dimension)         | Rxi-17 Sil MS 0.6 m x 0.25 mm i.d. x 0.25 µm (Restek, Bellefonte, PA, USA)  |
| Temperature Program                 | 50°C (0.50 min.), ramped 5°C/min. to 150°C (1.01 min.), ramped 2°C/min. to 200°C, ramped 50°C/min. to 300°C (15 min.)             |
| Modulation                          | 2 <sup>nd</sup> oven maintained +10°C relative to primary oven<br>4s with temperature maintained +15°C relative to secondary oven |
| Mass Spectrometer                   | LECO Pegasus® BT 4D   |
| Ion Source Temperature              | 250 °C  |
| Ionization Mode                     | El  |
| Mass Range (m/z)                    | 45-600  |
| Acquisition Rate                    | 15 spectra/s (1D); 200 spectra/s (2D)   |

## Why GCxGC for Urine Analysis?






Figure 1. A) GC-TOFMS TIC of Non-smoker's urine (NSU), B) extracted Ion Chromatogram (XIC) showing coeluting parabenic acid and D-(-)-citramalic acid, and C) corresponding Caliper, Peak True and Library Spectra.



## Confident Characterization: Similarity, Retention Index, and Mass Δ

### Spectral Similarity





Figure 4. NSU GCxGC-TOFMS Peak True spectra, library mass spectra, RI (Experimental and NIST calculated) and Mass Δ values for methylmalonic acid (A/B), and adipic acid (C/D).

| Name                                 | R.T. (s)       | Similarity | Name   | R.T. (s)        | Similarity |
|--------------------------------------|----------------|------------|--|-----------------|------------|
| Lactic acid, 2TMS                    | 588.007, 0.683 | 933        | Citric acid, 4TMS                            | 1968.312, 1.140 | 882        |
| Glycolic acid, 2TMS                  | 612.009, 0.736 | 915        | Methylcitric acid, 4TMS                      | 2004.12, 1.127  | 814        |
| Oxalic acid, 2TMS                    | 708.017, 0.894 | 927        | Adenine, 2TMS                                | 2020.12, 2.210  | 840        |
| 2-Methyl-3-hydroxybutyric acid, 2TMS | 816.025, 0.707 | 949        | m-Coumaric acid, 2TMS                        | 2028.12, 1.664  | 790        |
| 3-Hydroxyvaleric acid, 2TMS          | 836.027, 0.714 | 912        | 1,5-Anhydroxylitol, 4TMS                     | 2032.12, 1.066  | 826        |
| Guaiacol, TMS                        | 856.028, 1.085 | 805        | D-Fructose, MOX, 5TMS                        | 2120.13, 0.898  | 908        |
| 4-Hydroxybutanoic acid, 2TMS         | 876.03, 0.780  | 805        | L-Ascorbic acid, 2-O-methyl-3,5,6-tris-O-TMS | 2144.13, 1.464  | 794        |
| 2-Methylglutonic acid, TMS           | 884.031, 0.781 | 813        | d-Galactose, (1E)-MOX, 5TMS                  | 2164.13, 0.916  | 933        |
| Benzoic acid, TMS                    | 884.031, 1.104 | 932        | d-Galactose, (1Z)-MOX, 5TMS                  | 2216.14, 0.937  | 909        |
| Niacin, TMS                          | 960.037, 1.189 | 912        | d-Glucose, (1Z)-MOX, 5TMS                    | 2224.14, 1.064  | 844        |
| 1,2,3-Butanetriol, 3TMS              | 980.038, 0.634 | 935        | D-Mannitol, 6TMS                             | 2256.14, 0.829  | 915        |
| Glyceric acid, 3TMS                  | 1040.04, 0.729 | 920        | 1H-Indole-2-acetic acid, 2TMS                | 2256.14, 2.475  | 924        |
| Malonic acid, 3TMS                   | 1148.05, 0.798 | 902        | D-Sorbitol, 6TMS                             | 2276.14, 0.835  | 934        |
| 3-Aminoisobutyric acid, 3TMS         | 1228.06, 0.754 | 823        | Myo-Inositol, 6TMS                           | 2312.14, 0.849  | 923        |
| Amtranilic acid, TMS                 | 1256.06, 1.461 | 924        | D-Gluconic acid, 6TMS                        | 2432.15, 0.923  | 910        |
| Pyroglutamic acid, TMS               | 1272.06, 1.745 | 949        | Palmitic acid, TMS                           | 2448.16, 1.416  | 923        |
| Malic acid, 3TMS                     | 1284.06, 0.847 | 922        | Seyllo-Inositol, 6TMS                        | 2500.16, 0.903  | 938        |
| Uracil                               | 1320.07, 3.345 | 916        | Kynurenic acid, 2TMS                         | 2508.16, 2.538  | 852        |
| 3-Hydroxybenzoic acid, 2TMS          | 1392.07, 1.161 | 860        | N-Acetyl-D-glucosamine, MOX (anti), 4TMS     | 2604.17, 1.306  | 850        |
| Trigonelline TMS                     | 1420.07, 2.076 | 863        | N-Acetyl-D-glucosamine, MOX (syn), 4TMS      | 2624.17, 1.316  | 869        |
| 3-Hydroxyphenylacetic acid, 2TMS     | 1472.08, 1.215 | 937        | Stearic acid, TMS                            | 2808.18, 0.555  | 912        |
| 4-Hydroxybenzoic acid, 2TMS          | 1504.08, 1.247 | 877        | Xanthuric acid, 3TMS                         | 2824.19, 0.586  | 807        |
| 4-Hydroxybenzoic acid, 2TMS          | 1528.08, 1.360 | 904        | D-Lactose, MOX, 8TMS (isomer 2)              | 2956.2, 0.524   | 916        |
| Vanillyl alcohol, 2TMS               | 1532.08, 1.328 | 823        | Maltose, 8TMS (isomer 2)                     | 2964.2, 0.556   | 891        |
| Furoylglycine, TMS                   | 1556.08, 2.135 | 930        | D-(+)-Cellobiose, MOX, 8TMS (isomer 2)       | 2980.2, 0.581   | 873        |
| Vanillylmandelic acid, 3TMS          | 1644.09, 1.607 | 851        | Tryptophan, 4TMS                             | 2980.2, 1.315   | 805        |
| Levogluconic acid, 3TMS              | 1692.1, 1.174  | 919        | Maltose, 8TMS, isomer 1                      | 3302.2, 0.971   | 840        |
| Acetic acid, (E)-, 3TMS              | 1780.1, 1.276  | 831        | Sucrose, 8TMS                                | 3336.23, 0.977  | 846        |

## Mass Δ, Retention Index



## Smoker's Urine Results




Figure 5. Drugs in smoker's urine (SU).




Figure 6. Tobacco related compounds in SU.




Figure 7. SU GCxGC-TOFMS Peak True spectra, library mass spectra for cotinine (A/B), and trans-3'-hydroxycotinine (C/D).




Figure 8. Target Analyte Finding (TAF) processing method for rapid and robust identification of tobacco related compounds in comprehensive data files.

Table 5. Retention times, similarity values for drugs in SU

| Name                       | R.T. (s)       | Similarity |
|----------------------------|----------------|------------|
| Methyl salicylate, TMS     | 1132.05, 1.237 | 838        |
| Nudiflorine                | 1236.06, 3.209 | 811        |
| Bamethan, 2TMS             | 1308.06, 0.393 | 816        |
| Gabapentin lactam          | 1396.07, 2.796 | 938        |
| GabapentinLactam-TMS       | 1420.07, 1.613 | 944        |
| Ibuprofen, TMS             | 1492.08, 1.372 | 889        |
| Paracetamol, 2TMS          | 1520.08, 1.240 | 906        |
| Paracetamol, TMS           | 1748.1, 2.638  | 913        |
| Caffeine                   | 1944.12, 0.919 | 871        |
| Benadryl                   | 2060.12, 2.873 | 900        |
| Theophylline, TMS          | 2108.13, 3.314 | 793        |
| Methcathinone              | 2160.13, 0.336 | 814        |
| Naproxen, TMS              | 2524.16, 3.003 | 911        |
| Acetaminophen              | 3056.2, 0.477  | 892        |
| Pregabalin, O,N-N-tri-TMS- | 3084.21, 0.698 | 759        |

Table 6. Retention times, similarity values for tobacco related compounds in SU

| Name   | R.T. (s)       | Similarity |
|--|----------------|------------|
| 3-Pyridinol, TMS                               | 536.003, 0.864 | 812        |
| Benzene, 1,3-dichloro-                         | 536.003, 1.142 | 837        |
| Phenol, TMS                                    | 564.005, 0.857 | 898        |
| 2-Pyridinecarboxitrile                         | 584.007, 1.815 | 898        |
| Carbazole, 2,4,7-trimethyl-o-Cresol, TMS       | 672.014, 0.497 | 715        |
| o-Cresol, TMS                                  | 692.015, 0.886 | 780        |
| 3-Pyridinol, TMS                               | 704.016, 1.044 | 915        |
| p-Cresol, TMS derivative                       | 728.018, 0.932 | 886        |
| Naphthalene, 2,6-dimethyl-                     | 772.022, 0.675 | 752        |
| 3-Ethylphenol, TMS                             | 872.03, 0.957  | 781        |
| Catechol, 2TMS                                 | 1008.04, 0.881 | 923        |
| Pyrene, 1,9-dimethyl-                          | 1080.05, 0.849 | 724        |
| 4-Cyanophenol, TMS                             | 1104.05, 1.286 | 733        |
| 4-Methylcatechol, 2TMS                         | 1120.05, 0.892 | 908        |
| Hydroquinone, 2TMS                             | 1136.05, 0.899 | 786        |
| Cotinine                                       | 1648.09, 0.021 | 950        |
| Theobromine                                    | 1984.12, 1.660 | 918        |
| trans-3'-Hydroxycotinine, TMS                  | 2064.13, 3.100 | 928        |
| Theobromine, TMS derivative                    | 2224.14, 0.044 | 725        |
| 2-Hydroxy-3-methylanthraquinone, O-TMS         | 2408.15, 2.800 | 810        |
| 6-Hydroxy-α-methylnaphthaleneacetic acid, 2TMS | 2728.18, 1.694 | 897        |
| 4-Nitrophenyl-β-D-galacturonide, 3TMS          | 3028.2, 0.927  | 736        |

Table 7. TAF processing results for NSU and NSU. As expected increased quantities of tobacco related compounds were detected in smoker's urine

| Name   | R.T. (s)     | SU Area   | NSU Area  |
|--|--------------|-----------|-----------|
| 3-Pyridinol, TMS                               | 536 ± 0.868  | 371202153 | 310107452 |
| Benzene, 1,3-dichloro-                         | 536 ± 1.148  | 2650601   | 2571935   |
| Phenol, TMS                                    | 564 ± 0.860  | 277036716 | 149109915 |
| 2-Pyridinecarboxitrile                         | 584 ± 1.822  | 3093248   | 2150723   |
| Carbazole, 2,4,7-trimethyl-o-Cresol, TMS       | 676 ± 0.424  | 9471222   | 10467511  |
| o-Cresol, TMS                                  | 692 ± 0.905  | 1882395   | 494412    |
| p-Cresol, TMS                                  | 728 ± 0.948  | 401969495 | 429516716 |
| Naphthalene, 2,6-dimethyl-                     | 772 ± 0.702  | 560201    | 2488416   |
| 3-Ethylphenol, TMS                             | 872 ± 0.985  | 1525954   | 1005120   |
| Catechol, 2TMS                                 | 1008 ± 0.921 | 523485419 | 416668899 |
| Pyrene, 1,9-dimethyl-                          | 1080 ± 0.896 | 1643952   | 890189    |
| 4-Cyanophenol, TMS                             | 1104 ± 1.498 | 19191     | 15717     |
| 4-Methylcatechol, 2TMS                         | 1120 ± 0.940 | 104992508 | 82018098  |
| Hydroquinone, 2TMS                             | 1136 ± 0.936 | 16412520  | 9694548   |
| Cotinine                                       | 1648 ± 0.122 | 20447248  | 430248    |
| Theobromine                                    | 1984 ± 1.789 | 9072840   | 17023415  |
| trans-3'-Hydroxycotinine, TMS                  | 2060 ± 3.020 | 13862     | 8859      |
| 2-Hydroxy-3-methylanthraquinone, TMS           | 2404 ± 2.778 | 794       | 4223      |
| 6-Hydroxy-α-methylnaphthaleneacetic acid, 2TMS | 2728 ± 1.873 | 6364834   | 14267     |
| 4-Nitrophenyl-β-D-galacturonide, 3TMS          | 3032 ± 1.028 | 79003259  | 18560460  |

## Summary

- The Pegasus BT 4D facilitates fast and confident compound identification through enhanced two-dimensional chromatographic resolution and high performance TOFMS.
- Robust compound identification was obtained through spectral similarity searches of large, well-established databases, mass Δ calculations and retention index filtering.
- Comprehensive data can be processed via non-targeted Peak Find or rapid Target Analyte Finding.



%Error = 0.54%