

Spice Wars—Are You Battle Ready? Analysis of Synthetic Cannabinoids via Gas Chromatography–High Resolution Time-of-Flight Mass Spectrometry

David E. Alonso¹, John Rorabeck², and Joe Binkley¹ | ¹LECO Corporation, St. Joseph, MI; ²Berrien County Forensic Laboratory, Berrien Springs, MI

Background

- Since the mid-2000s, synthetic drugs have been at the forefront of a world-wide market in “legal high” mind-altering substances sold to customers without proper manufacturing protocol, quality controls, general safety studies or dosing information.
- They are available over-the-counter or via the internet as plant fertilizer, incense, potpourri, or bath salts.
- They are not typically ordered on a routine lab test panel.
- They include the following series of compounds: JWH, CP, HU, AM, WIN, XLR, and UR (Fig. 1).



Fig. 1: Generations of Synthetic Cannabinoids

Mr. Nice Guy

- Synthetic cannabinoids, such as Mr. Nice Guy (JWH-018 & JWH-073), burst onto the world scene as first-generation synthetic cannabinoids. Figure 4 shows the analytical ion chromatogram (AIC) for Mr. Nice Guy, as well as the Peak True (Deconvoluted) mass spectra for JWH-073 & JWH-018.




Fig. 4: A) Analytical Ion Chromatogram (AIC) of Mr. Nice Guy Extract. Peak True (Deconvoluted) Mass Spectra for B) JWH-073 and C) JWH-018

- In 2013, a seized drug packet was found to contain XLR-11, XLR-11 Isomer, and an unknown (Fig. 5). Acquisition of EI and CI-HRT data and a subsequent database search resulted in identification of the unknown as AB-Pinaca (Fig. 6).

Case 1: Botanical (Pokeweed)

