

Untargeted Investigation of Non-Alcoholic Fatty Liver Disease Using Effective Mutiplatform GC-MS Instrumentation

David E. Alonso¹, Xiang Zhang², Todd Richards¹, and Joe Binkley¹ | ¹LECO Corporation, St. Joseph, Michigan | ²Department of Chemistry and Biochemistry, University of Louisville, Louisville, Kentucky

Introduction

Background

- Non-alcoholic fatty liver disease (NAFLD) is a major medical issue in developed countries.
- NAFLD is the most common liver disease and parallels the obesity epidemic in the US.
- The purpose of this investigation is to utilize high performance Gas Chromatography Time-of-Flight Mass Spectrometry (GC-TOFMS) for increased analysis throughput of complex samples, followed by comprehensive, multidimensional chromatography and high resolution Time-of-Flight mass spectrometry (GCxGC-HRT) for discovery and confirmation of NAFLD biomarkers.

Methods

Samples

Rodents were fed different levels of copper with either distilled water or 30% fructose (w/v) for a period of several weeks. The animals were then euthanized and liver/plasma samples collected for TOFMS analysis. Metabolites were extracted with a water/methanol solvent mixture, dried extensively, and derivatized via a standard two-step procedure (MEOX, MTBSTFA). This methodology minimized sugar interferences (e.g., monosaccharides and disaccharides), and produced stable derivatives of alcohols, lactones, acids, amino acids, diacids, fatty acids, and sterols that could be easily analyzed by GC-TOFMS.

Instrumentation

Samples were analyzed using a combination of nominal and high resolution GC-TOFMS and GCxGC-HRT. Complementary electron and chemical ionization methods were also employed to assist in identification of metabolites. Acquired data were processed using untargeted methods for peak determination, spectral deconvolution, accurate mass formula determination, and spectral similarity searches for confident metabolite identifications.

Table 1. GC-TOFMS Parameters

Gas Chromatograph	Agilent 7890 with LECO L-Pal3 Autosampler
Injection	1 µL, Split 20:1; 270 °C
Carrier Gas	He @ 0.8 mL/min, Constant Flow
Column	Rxi-5 Sil MS, 20 m x 0.18 mm i.d. x 0.18 µm (Restek, Bellefonte, PA, USA)
Temperature Program	60 °C (0.50 min), ramped 36 °C/min to 320 °C (3 min)
Mass Spectrometer	LECO Pegasus BT
Ion Source Temperature	250 °C
Ionization Mode	EI
Mass Range (m/z)	45-600
Acquisition Rate	10 spectra/s



Figure 1. Pegasus® BT GC-TOFMS.

Table 2. GC and GCxGC-HRT Parameters

Gas Chromatograph	Agilent 7890 with LECO L-Pal3 Autosampler
Injection	1 µL, Split 20:1; 280 °C (1 µL, Splitless for CI)
Carrier Gas	He @ 1.0 mL/min, Constant Flow
Column 1	Rxi-5 MS, 30 m x 0.25 mm i.d. x 0.25 µm (Restek, Bellefonte, PA, USA)
Column 2	Rxi-17 Sil MS, 0.60 m x 0.25 mm i.d. x 0.25 µm (Restek, Bellefonte, PA, USA)
Temperature Program	60 °C (0.50 min), ramped 5 °C/min to 270 °C (6 min) Secondary oven maintained +10 °C relative to primary oven
Thermal Modulation (GCxGC)	4s with temperature maintained +15 °C relative to secondary oven
Mass Spectrometer	LECO Pegasus HRT
Ion Source Temperature	250 °C (EI); 200 °C (CI)
Ionization Mode	EI and CI (Reagent Gas: 5% NH ₃ in CH ₄)
Mass Range (m/z)	30-510 (EI); 60-1500 (CI)
Acquisition Rate	10 spectra/s (200 spectra/s GCxGC)



Figure 2. Pegasus GC-HRT 4D.

Throughput: GC-TOFMS




Figure 3. A) eXtracted Ion Chromatogram (XIC) and B) TIC of a liver sample.

Figure 4. A) ChromaTOF® brand software deconvolution results - liver sample. Peak True (deconvoluted) and library spectra for aconitic acid (B,C) and lysine (D,E).

Table 3. Representative analytes in liver sample (Similarity Ave. = 886)

Metabolite	Formula	R.T.(Sec)	Similarity
Lactic acid, 2TBDMS derivative	C ₁₁ H ₂₀ O ₅ Si ₂	273	866
Glycolic acid, 2TBDMS derivative	C ₁₀ H ₁₈ O ₅ Si ₂	276	865
Alanine dTBDMS	C ₁₁ H ₂₀ NO ₂ Si ₂	283	894
Glycine, 2TBDMS derivative	C ₁₀ H ₁₈ NO ₂ Si ₂	287	909
3-Hydroxybutyric acid, (R), 2TBDMS derivative	C ₁₁ H ₂₀ O ₅ Si ₂	294	911
L-2-Aminobutyric acid, 2TBDMS derivative	C ₁₀ H ₁₈ NO ₂ Si ₂	297	889
D-Alanine, 2TBDMS derivative	C ₁₁ H ₂₀ NO ₂ Si ₂	303	856
Dimethylglycine, 2TBDMS derivative	C ₁₁ H ₂₀ NO ₂ Si ₂	306	900
Valine dTBDMS	C ₁₁ H ₂₀ NO ₂ Si ₂	307	864
L-Leucine, 2TBDMS derivative	C ₁₃ H ₂₄ NO ₂ Si ₂	314	937
Niacinamide, 2TBDMS derivative	C ₁₁ H ₁₆ NO ₂ Si ₂	318	922
Isoleucine, 2TBDMS derivative	C ₁₃ H ₂₄ NO ₂ Si ₂	320	943
4-Aminobutyric acid, 2TBDMS derivative	C ₁₀ H ₁₈ NO ₂ Si ₂	325	831
Uracil, 2TBDMS derivative	C ₁₁ H ₁₆ NO ₂ Si ₂	327	929
L-Proline, 2TBDMS derivative	C ₁₁ H ₂₀ NO ₂ Si ₂	328	918
Racemic acid, 2TBDMS derivative	C ₁₁ H ₂₀ O ₅ Si ₂	331	864
D-Pyroglytamamic acid, 2TBDMS derivative	C ₁₁ H ₂₀ NO ₂ Si ₂	360	854
Mel, (2TBDMS)	C ₁₁ H ₁₆ NO ₂ Si ₂	363	868

Statistical Analysis: Features → Metabolites




Figure 5. XCMS online Cloud Plot (m/z vs retention time) showing a metabolite feature with a retention time of 7.55 minutes and EIC m/z = 440.39.

Figure 6. ChromaTOF - A) XIC 440.40; B) generated formulas for fragment and molecular ions; C, D) Peak True and library match data for histidine.

Table 4. Partial list of potential NAFLD biomarkers identified using XCMS online and ChromaTOF. Metabolites A-D were upregulated in liver samples, while citric acid was found in greater quantities in plasma. Stearic Acid was tentatively identified with the marginally acceptable similarity score of 703/1000

Feature	m/z	BT Similarity	Metabolite
A	341.28	703	Stearic Acid?
B	418.36	897	Aspartic Acid
C	431.43	917	Lysine
D	440.39	895	Histidine
E	459.33	726	Citric Acid




Figure 7. XCMS online Results Table showing a feature with a retention time of 6.52 minutes and EIC m/z = 418.36. ChromaTOF processing resulted in the identification of the metabolite aspartic acid.

Robust Characterization and Confirmation: GC-HRT




Figure 8. Peak True EI-HRT, library match, and CI-HRT data for stearic acid.

Discovery: GCxGC-HRT

- Enhanced Chromatographic and Mass Spectral Resolution
- Group Clustering – Structured Chromatograms
- Removal of background interferences (e.g., column bleed, solvents, etc.)
- Improved Characterization of Compounds in Complex Matrices





Figure 9. A) GC-HRT and B) GC-HRT 4D data. A table illustrating the benefits of combining enhanced chromatographic and mass spectral resolution for confident identification of metabolites differentiating the different groups (e.g., fumaric acid).




Figure 10. A) GC-HRT and B) GC-HRT 4D data for pantothenic acid (Vitamin B₅).

Summary

- High performance GC-TOFMS was used for increased analysis throughput of complex samples.
- GC-HRT facilitated confident identification of metabolites using complementary high resolution mass spectrometry ionization methods for formula determinations for fragment, molecular, and adduct ions.
- The combination of enhanced chromatographic and mass spectral resolution (Figure 11) led to superior analysis of complex matrices by transforming unknowns to known compounds in metabolomics data. Compounds were confidently identified as potential biomarkers for NAFLD as a direct result of this methodology.




Figure 11. GC-HRT 4D data for citric acid.