# Utilization of GCxGC-TOFMS to Screen for Potential Metabolite Differences in Pooled Samples from Lean, Fat, and Obese Zucker Rat Plasma John Heim, Jeff Patrick, and Joe Binkley • Saint Joseph, MI USA

### **OBJECTIVES**

Demonstrate the enhanced detection and increased peak capacity available with GCxGC-TOFMS for complex metabolite comparisons.

Show how comprehensive two-dimensional chromatography resolves analytes that would wise totally coelute and be buried by high concentration analytes in one-dimensional otherwise totally

Illustrate the value of overlaid pooled sample chromatographic and TOFMS data comparison to discover distinct analyte variations

Show how GCxGC-TOFMS can be used as a fast screening tool of large sample pools to identify potential biomarkers in lean, fat, and diabetic obese Zucker rats.

### INTRODUCTION

This GCxGC-TOFMS research presents the comparison of pooled derivatized plasma samples from lean, Zucker rats. The small molecule metabolite profiles from the different sample pools were compared for similarity in mass spectral identification, retention time, and analyte variation by peak area. The enhanced chromatographic resolution and peak capacity available with comprehensive two-dimensional chromatography provided by the fast acquisition of TOFMS, up to 500Hz, that will successfully acquire the data density needed to fully characterize low levels of metabolites common in complex sample matrices such as blood plasma. These complex and data rich files were processed with deconvolution algorithms which deliver qualitative identification.

GCxGC-TOFMS was used to evaluate pooled trimethylsilyl-derivatized rat plasma samples in various states of health. Plasma from each of the strains of rat (Zucker Lean, Zucker Fatty and Zucker Diabetic Obese from Charles River Labs) was derivatized with BSTFA after the proteins were removed. Sample analysis was conducted on the lean, fat, and obese pooled samples utilizing a GCxGC separation followed by TOFMS detection at an acquisition rate of 150 spectra per second. The raw data was processed and mass spectral searches conducted using Max Planck, Fiehn Rtx5, and NIST libraries. The data mining strategy compared the combined peak tables, mass spectral similarities, and overlaid chromatograms to discover significant metabolite variations between the lean, fat, and obese sample pools.

Rat Plasma was obtained from Bioreclamation. The GCxGC-TOFMS analysis was performed on derivatized pooled samples from various states of health in lean, fat, and diabetic obese Zucker rats. This exploratory research was conducted as a screening tool to characterize and measure semi-quantitative variations in metabolites between the sample pools. The raw data files were first processed then compared by the chromatographic first and second dimension retention time, and mass spectral similarity. A total of twelve metabolites between the different sample pools were found to have significant variation measured by their chromatographic peak areas. This exploratory research demonstrates the favorable and practical applicability of GCxGC-TOFMS as a simple and time saving screening tool to discover potential metabolic biomarkers which can be utilized to distinguish between healthy and disease state sample pools.

### **EXPERIMENTAL METHODS**

#### **Sample Preparation**

- Proteins were removed from lean, fat and obese Zucker rat plasma with a 5000 MWCO (molecular weight cut-off) filter.
- 100µL samples were dried in a Savant SpeedVac.
- The dried residue was suspended in 100µL BSTFA (N,O-bis(trimethylsilyl)trifluoroacetamide) and incubated overnight at 60°C.
- GCxGC-TOFMS analysis was performed on the pooled samples from various states of health (lean, fat and diabetic obese).

### GCxGC-TOFMS Analysis Parameters

- Gas Chromatograph: Agilent 7890 equipped with a LECO dual stage guad jet thermal modulator
- Primary Column: 30 m x 0.25 mm id. x 0.25 µm film thickness Rxi-5Sil MS (Restek Corp., Bellefonte, PA) Secondary Column: 1.24 m x 0.15 mm id. x 0.15 µm film thickness Rxi-17Sil MS (Restek Corp., Bellefonte,
- PA)
- Carrier gas: Helium set @ 1.5 mL/min
- Injection: 1 µL splitless @ 260°C
- Primary column temperature program: Initial temperature 70°C for 0.5min., then ramped 6°C/min. to 305°C for 5 min.
- Secondary column temperature program: Parallel ramp offset by +5°C
- Total run time: 44.67 min.

#### Mass Spectrometer: LECO Pegasus<sup>®</sup> 4D

- Mass Range: 45 750 m/z
- Acquisition rate: 150 spectra/s
- Ion source temperature: 240°C
- Detector Voltage: 1650 V
- Electron Energy: -70 eV





**Figure 1.** The figure above shows the overlaid linear unique mass ion chromatogram of Leucine, N,O-TMS for the lean, fat, and diabetic obese pooled sample GCxGC-TOFMS analysis. The mass spectra show that leucine was not found in the lean pooled sample while the mass spectra for the fat, and obese pooled sample are identified as leucine. The 3D surface plot chromatograms for lean, fat, and obese pooled samples indicate that leucine is not found in the lean pooled sample and is upregulated in the fat and diabetic obese sample pools.







Figure 2. Figure 2 above shows the overlaid linear unique mass ion chromatogram of TMS-ester-oxobutyl-glycine for the lean, fat and diabetic obese pooled sample GCxGC-TOFMS analysis. The mass spectra show that TMS-ester-oxobutyle glycine was not found in the lean pooled sample. However, it is identified in both the fat and obese pools. The 3D surface plots of the lean, fat, and obese pooled samples indicate that TMS-oxobutyl-glycine is upregulated in the fat and diabetic obese sample pools and not present in the lean sample pool.



**Delivering the Right Results** 

### **GCxGC-TOFMS POOLED SAMPLE METABOLITE COMPARISON RESULTS**





**Figure 3.** A two dimensional contour plot of the diabetic obese pooled sample is shown above in figure 3. A total of 438 peaks were found with a S/N ratio of greater than or equal to 100. The GCxGC contour plot visually illustrates the enhanced peak capacity and resolution that is not possible by a one dimension separation.







in all 3 sample pools. The surface plots for lean, fat, and obese pooled samples indicate that N,O-TMS-Proline is upregulated from the lean to the fat sample pool and then downregulated in the diabetic obese sample pool.

990) 1

**Figure 5.** Figure 5 shows the overlaid linear unique mass ion chromatogram of Penta-TMS-Glucopyranose for the lean, fat, and diabetic obese pooled sample GCxGC-TOFMS analysis. The mass spectra show that Penta-TMS-Glucopyranose was found in all 3 sample pools. The surface plots for lean, fat, and obese pooled samples indicate the expected trend that Penta-TMS-Glucopyranose is upregulated from the lean to the fat and the obese diabetic sample pools.



Figure 6. Figure 6 shows the overlaid linear unique mass ion chromatogram of Tris-TMS-L-Threonine for the lean, fat and diabetic obese pooled sample GCxGC-TOFMS analysis. The mass spectra show that Tris-TMS-L-Threonine was found in all 3 sample pools. The surface plots for the lean, fat, and obese pooled samples indicate an interesting trend that Tris-TMS-L-Threonine is upregulated from the lean to the fat sample pool and downregulated to the least concentration in the diabetic obese sample pool.

### **GCxGC ENHANCED PEAK CAPACITY**



Figure 7. The three dimensional surface plot chromatogram shown above illustrates the increased peak capacity and resolution of comprehensive two-dimensional chromatography. Notice the separation of trace level peaks that are found in the second dimension plane. The white total ion chromatogram of the 1<sup>st</sup> dimension separation is shown in the top right-hand section of figure 7. This example clearly illustrates that GCxGC resolves and separates components that would otherwise be coeluted and buried under high concentration analytes by a one dimensional separation.

## CONCLUSIONS

This study highlights the use of GCxGC-TOFMS analysis of pooled Zucker rat plasma samples from different states of health. The chromatographic and time-of-flight mass spectral comparisons between sample pools establishes this analytical approach as a valuable screening tool for discovery of metabolite variation across different classes. The analysis of derivatized complex blood plasma demonstrates the enhanced detection and increased peak capacity available with GCxGC-TOFMS using the LECO Pegasus 4D instrument. The GCxGC-TOFMS results for the lean, fat, and obese sample pools demonstrates how GCxGC separates and resolves components that would otherwise totally coelute and be buried under high concentration analytes in one-dimensional chromatography. TOFMS provides non-skewed mass spectral data and the fast acquisition necessary to facilitate deconvolution algorithms that yield high quality library searchable mass spectra. This proof of concept research shows a fast, efficient, and simple comparative analysis method for screening large sample pools to distinguish metabolite variations as potential biomarkers for healthy and disease state populations. GCxGC-TOFMS is a valuable time and resource savings instrumental option for the comparative analysis of pooled samples from varying states of health to quickly and efficiently determine significant metabolic differences

For further information regarding this study contact the authors at john\_heim@leco.com