

# Smoker's and Non-Smoker's Urine Comparison Using Comprehensive Two-Dimensional Gas Chromatography High Performance Time-of-Flight Mass Spectrometry

David E. Alonso\* and Joseph E. Binkley | LECO Corporation, St. Joseph, Michigan U.S.A.

## Introduction

Urine is a favored biofluid for diagnostic testing (Urinalysis) because it is non-invasive and large volumes are easily obtained. In addition, urine is relatively free from interfering proteins and lipids, and it tends to "hold" high concentrations of drugs and metabolites over extended periods of time. Modern, routine clinical tests include the determination of specific gravity, measurement of glucose, nitrates, etc.

In this study, a novel analytical approach was utilized for the effective characterization of compounds in two standard reference materials, NIST smoker's and non-smoker's urine.

## Objectives

- Implement the use of enhanced, comprehensive two-dimensional gas chromatography (GCxGC) for the separation of compounds in urine
- Use a benchtop, high performance time-of-flight mass spectrometer and powerful processing software to quickly and confidently identify compounds in urine
- Use software tools to compare smoker's and non-smoker's urine

## Sample Preparation

- 600  $\mu$ L of urine treated with urease ( $37^{\circ}\text{C}$ , 15 min)
- The mixture was vortexed (2 min) and then centrifuged (12,000 g for 10 min)
- 200  $\mu$ L of supernatant was transferred to a 2mL GC vial and evaporated to dryness (Speed Vac)
- The dry material was derivatized using a two-step procedure
  - Methoximation (20  $\mu$ L of MEOX,  $80^{\circ}\text{C}$ , 30 min)
  - Silylation (75  $\mu$ L MSTFA,  $80^{\circ}\text{C}$ , 30 min)

## Technology




Table 1. Instrument acquisition parameters

Gas Chromatograph	Agilent 7890, LECO Dual Stage Quad Jet Modulator & L-PAL 3 Autosampler
Injection	1 $\mu$ L, Split 20:1; $280^{\circ}\text{C}$
Carrier Gas	He @ 1.4 ml/min, Constant Flow
Columns (1 <sup>st</sup> Dimension)	Rxi-5 MS, 30 m x 0.25 mm i.d. x 0.25 $\mu\text{m}$ (Restek, Bellefonte, PA, USA)
(2 <sup>nd</sup> Dimension)	Rxi-17 Sil MS 0.6 m x 0.25 mm i.d. x 0.25 $\mu\text{m}$ (Restek, Bellefonte, PA, USA)
Temperature Program	50 °C [0.50 min], ramped 5 °C/min to 150 °C (1.01 min), ramped 2 °C/min to 200 °C, ramped 50 °C/min to 300 °C (15 min)
Modulation	2 <sup>nd</sup> oven maintained +10 °C relative to primary oven
Modulation	4s with temperature maintained +15 °C relative to secondary oven
Mass Spectrometer	LECO Pegasus BT 4D
Ion Source Temperature	250 °C
Ionization Mode	El
Mass Range (m/z)	45-600
Acquisition Rate	15 spectra/s (1D); 200 spectra/s (2D)

## Why GCxGC for Urine Analysis?






Table 2. Comparison of GC and GCxGC-TOFMS spectral similarity values for some acids in NSU (Unknowns → Knowns)

GC-TOFMS			
Name	R.T. (s)	Similarity	Mass Δ (Da)
D-(+)-Citramalic acid, 3TMS	1255.77	756	
Parabanic acid, 2TMS		Not Found	
Kojic acid, 2TMS	1264.51	570	
Quinolinic acid, 2TMS	1721.83	326	
Orotic Acid, 3TMS	1778.33	660	
Homovanillic Acid, 2TMS	1814.9	449	
Hippuric acid, TMS	1847.01	649	
Vanillylmandelic acid, 3TMS	2102.38	681	
Panthothenic acid, 3TMS	2369.86	537	
Caffeic acid, 3TMS	2702.75	533	

Name	R.T. (s)	Similarity	Mass Δ (Da)
D-(+)-Citramalic acid, 3TMS	1256.06, 0.802	864	N/A
Parabanic acid, 2TMS	1256.06, 1.199	855	0.01
Kojic acid, 2TMS	1264.06, 1.020	818	-0.01
Quinolinic acid, 2TMS	1720.1, 2.088	941	0.01
Orotic Acid, 3TMS	1776.1, 1.357	807	0.02
Homovanillic Acid, 2TMS	1814.9, 449	895	0.01
Hippuric acid, TMS	1844.12, 2.811	953	0.02
Vanillylmandelic acid, 3TMS	2096.13, 1.537	860	0.02
Panthothenic acid, 3TMS	2364.15, 1.395	912	N/A
Caffeic acid, 3TMS	2696.18, 1.615	865	0.02

## Confident Characterization: Similarity, Retention Index, and Mass Δ

### Spectral Similarity




Table 3. Representative list of compounds in NSU with retention times and spectral similarity values

Name	R.T. (s)	Similarity	Name	R.T. (s)	Similarity
Lactic Acid, 2TMS	588.007, 0.683	933	Citric acid, 4TMS	1968.12, 1.140	882
Glycolic acid, 2TMS	612.009, 0.736	915	Methyltric acid, 4TMS	2004.12, 1.127	814
Oxalic acid, 2TMS	708.017, 0.893	927	Adenine, 2TMS	2020.12, 2.210	840
2-Methyl-3-hydroxybutyric acid, 2TMS	816.025, 0.707	949	1,5-Anhydroyhexitol, 4TMS	2032.12, 1.066	826
3-Hydroxyvaleric acid, 2TMS	856.028, 1.085	805	D-Fructose, MOX, 5TMS	2120.13, 0.898	908
4-Hydroxybutanoic acid, 2TMS	876.030, 0.780	805	L-Ascorbic acid, 2-O-methyl-3,5,6-tris(2,4-dimethylphenyl)-Abundance	2144.13, 1.464	794
2-Methyloctanoic acid, 2TMS	884.031, 0.783	813	d-Galactose, (1E)-MOX, 5TMS	2164.13, 0.916	933
Benzoic Acid, TMS	890.037, 0.731	932	d-Galactose, (1Z)-MOX, 5TMS	2216.14, 0.937	909
Nicotinic Acid, 2TMS	900.037, 0.734	912	D-Mannose, 2TMS	2224.14, 1.044	944
1,2,3,4-Tetrahydro-2H-pyran, 2TMS	909.037, 0.734	915	1-H-indole-3-acetic acid, 2TMS	2256.14, 2.475	934
Glycric acid, 3TMS	1040.04, 0.729	920	D-Sorbitol, 6TMS	2276.14, 0.835	934
Malonic acid, 3TMS	1148.05, 0.798	902	Myo-Inositol, 6TMS	2312.14, 0.849	923
3-Aminobutyric acid, 3TMS	1228.06, 0.754	823	Anthrancic acid, TMS	2423.06, 0.923	910
Euroglutamic acid, TMS	1256.06, 1.461	924	Pyroglutamic acid, TMS	2472.06, 1.745	949
Malic acid, 3TMS	1284.06, 0.847	922	Palmitic Acid, TMS	2448.16, 1.416	923
3-Hydroxybenzoic acid, 2TMS	1392.07, 1.166	860	N-Acetyl-D-glucosamine, MOX (anti), 4TMS	2508.16, 2.538	852
Trigonelline, TMS	1472.07, 0.815	863	N-Acetyl-D-glucosamine, MOX (syn), 4TMS	2624.17, 1.316	869
3-Hydroxypropanoic acid, 2TMS	1472.08, 1.315	937	Stearic acid, 2TMS	2688.16, 0.935	912
4-Hydroxybutyric acid, 2TMS	1528.08, 1.360	904	D-Lactose, MOX, 8TMS (isomer 2)	2956.05, 0.526	916
Vanillyl alcohol, 2TMS	1528.08, 1.328	823	Malto-8TMS (isomer 2)	2964.2, 0.556	891
Euroglutone, TMS	1556.08, 2.135	930	D(+)-Cellobiose, MOX, 8TMS (isomer 2)	2980.2, 0.581	873
Vanillylmalonic acid, 3TMS	1644.09, 1.607	851	Tryptophan, 4TMS	3000.2, 2.315	805
Levoglucosan, 3TMS	1692.1, 1.174	919	Malto-8TMS, isomer 1	3304.22, 0.971	840
Aconitic acid, (E), 3TMS	1780.11, 1.276	831	Sucrose, 8TMS	3336.23, 0.977	846

X %Error = 0.54

Table 4. Comparison of experimental and NIST RI values for diacids in urine

Name	R.T. (s)	Similarity	Mass Δ (Da)	Exp. RI	NIST RI
Oxalic acid, 2TMS	708.017, 0.893	927	N/A	1142	1136
Methylmalonic acid, 2TMS	844.026, 0.825	778	0.02	1224	1223
Succinic acid, 2TMS	1004.04, 0.896	885	0.01	1323	1321
Methylsuccinic acid, 2TMS	1024.04, 0.861	907	N/A	1336	1331
Fumaric acid, 2TMS	1052.04, 0.823	915	N/A	1354	1353
Iaconic acid, 2TMS	1052.04, 0.942	816	N/A	1354	1398
Methylmaleic acid, 2TMS	1064.05, 0.972	876	N/A	1362	1386
3-Methylmalonic acid, 2TMS	1172.05, 0.892	906	N/A	1432	
Adipic acid, 2TMS	1296.06, 0.993	854	-0.02	1514	1514
3-Methyladipic acid, 2TMS	1344.07, 1.021	849	N/A		