Development of Comprehensive Steroid Analysis Methods by GCxGC-HR-TOFMS

<u>Viatcheslav Artaev¹</u>, Jonelle Shiel¹, David Alonso¹, Georgy Tikhonov¹, Albert Lebedev² | ¹LECO Corporation, St. Joseph, MI, USA ² M.V. Lomonosov Moscow State University, Chemistry Department, Moscow, Russian Federation

Introduction

Identification and quantitation of steroid metabolites in biological samples are essential for screening for various hormonal disorders. Many of these metabolites are closely related isomers and cannot be easily separated by LC or one-dimensional GC due to their chemical and structural similarity, which complicates or even prevents reliable analyte assignment^[1]. Comprehensive two-dimensional gas chromatography (GCxGC) coupled to a high resolution, high mass accuracy time-of-flight mass spectrometer (HR-TOFMS) provides dramatically enhanced chromatographic separation, increase in sensitivity, and reliable detection and identification of analytes of interest.

Thirty-three steroids from different classes (progestogens, androgens, estrogens, glucocorticoids, mineralocorticoids) were derivatized and analyzed using GCxGC coupled to HR-TOFMS, to achieve detection and reliable identification in complex matrices.

Methods

The analytes were run as individual components first, and the resulting high resolution mass spectra were used to create a custom accurate mass library (AML). The AML assisted in the identification of steroids in the biological matrices. In addition, the standard mixture of the steroid was used for development of a chromatographic method for achieving the most efficient multidimensional separation by applying Simply GCxGC[®] software tool.^[2] ChromaTOF[®] brand software (LECO, St. Joseph, MI) was used for instrument control, data acquisition, AML creation, peak finding, and compound identification.



Results



Steroid Standards Sample Preparation

For this study we acquired 32 steroid standards from Steraloids and one standard from Fisher Scientific.

Standard Preparation:

Prepare a 2 mg/ml solution in Methanol; Pipette 10 μ l into an autosampler vial; Speed Vac to dryness for ~ 30 minutes

Derivatization:

• Add 20 μ l Methoxyamine HCl in pyridine (20 mg/ml); Heat and agitate at 80 °C for 1 hour; Add 80 μ I MSTFA +1% TMCS (purchased from Thermo Scientific); Heat and agitate at 100 °C for 1 hour







Figure 4. 2D Chromatogram Contour Plot of the derivatized mixture.



Abbreviation	Trivial Name	Abbreviation	Trivial Name
SS (ISTD)	Stigmasterol	MP (ISTD)	Medroxyprogesterone
ANDRO	Androsterone	THA	Tetrahydro-11-dehydrocorticosterone
ETIO	Etiocholanolone	5a-THB	5a-Tetrahydrocorticosterone
DHA	Dehydroepiandrosterone	THF	Tetrahydrocortisol
11-OXO-ETIO	11-oxo-Etiocholanolone	5a-THF	5a-Tetrahydrocortisol
17β-Estradiol	17β-Estradiol	a-Cortolone	a-Cortolone
17-HP	17-Hydroxypregnanolone	β-Cortol	β-Cortol
11β-OH-ANDRO	11β-Hydroxyandrosterone	β-Cortolone	β-Cortolone
16a-OH-DHA	16a-Hydroxy-DHEA	Cortisone	Cortisone
PT	Pregnanetriol	Cortisol	Cortisol
5-AT	Androstentriol	20β-DHE	20β-Dihydrocortisone
THS	Tetrahydro-11-deoxycortisol	20a-DHE	20a-Dihydrocortisone
THDOC	Tetrahydrodeoxycorticosterone	20β-DHF	20β-Dihydrocortisol
Estriol	Estriol	6β-OH-F	6β-Hydroxycortisol
PT'ONE	Pregnanetriolone	18-OH-F	18-Hydroxycortisol
5-PT	Pregnentriol, 5-PT	20a-DHF	20a-Dihydrocortisol
THE	Tetrahydrocortisone		



Figure 5. Contour 2D Chromatogram Plot region demonstrating an example of separation of the isomeric compounds.





Figure 1. Pegasus[®] HRT⁺ 4D (LECO Corp., St. Joseph, MI) – GCxGC-HR-TOFMS system used in this study.

Gas Chromatograph - Agilent 7890				
Injection	1μL, Split 100:1, 250 °C			
Carrier Gas	He, 1.4 mL/min			
Temperature	200 °C (0.5 min) – 300 °C at 5 °C/min			
GCxGC – LECO Cryoge	enic Thermal Modulator			
Columns	1D: 15 m x 0.250 mm x 0.25 μm HP-1MS			
	2D: 2 m x 0.250 mm x 0.25 μm BPX-50			
	2D: 1.75 m in GC Oven, 0.10 m in Modulator, 0.15 m in 2D Oven			
	Guard Column 1.4 m x 0.250 mm x 0 μ m Uncoated			
	Guard Column: 0.80 m in 2D Oven, 0.60 m in Transfer Line			
Temperaure	2D Oven: +13 °C, Modulator: +15 °C			
Modulation Period	3 seconds, Hot Pulse: 0.9 second			
Mass Spectrometer –	LECO Pegasus HRT+ 4D (R=25K @ 219 M/Z)			
Transferline Temperature	300 °C			
Ion Source Temperature	250 °C			
Spectra Acquisition Rate	200 spectra/second			
Mass Range, m/z	40-1000			

Summary

- AML library was created by running individual derivatized steroid samples.
- GCxGC Method was developed allowing reliable separation of all 33 steroids in the sample.
- All 33 analytes were found in the mixture and positively assigned using accurate mass confirmation for molecular ions (when available) and major fragment ions, RI, and AML/NIST library matches.

Mixture Preparation

Mixture Preparation:

Pipette 10 μ l of the 2 mg/ml of each of the standards into an autosampler vial; Speed Vac to dryness for ~ 1 hour

Derivatization MSTFA:

• Add 60 μ l Methoxyamine HCl in pyridine (20 mg/ml); Heat and agitate at 80 °C for 1 hour; Add 300 μ l MSTFA +1% TMCS (purchased from Thermo); Heat and agitate at 100 °C for 1 hour

Derivatization TMSI:

• Add 60 μ l Methoxyamine HCl in pyridine (20 mg/ml); Heat and agitate at 80 °C for 1 hour; Speed Vac to dryness for ~ 1 hour; Add 350 μ l TMSI (purchased from Thermo Scientific); Heat and agitate at 100 °C for 16 hours



Future Work

- Optimize derivatization procedures of the steroids in the matrix (urine)
- Validate the GCxGC-HR-TOFMS analysis method by using urine samples spiked with the standards mixture

References

^[1] Comprehensive Steroid Analysis By GCxGC-TOFMS, Michael Groessl et al, 288477, ThOG, Proceedings of the 65th ASMS **Conference on Mass Spectrometry and Allied Topics,** Indianapolis, IN, June 4-8, 2017.

^[2] https://www.leco.com/simply-gcxgc

