# Fragmentation Trees for Automated de novo Interpretation of Impure Electron Ionization Spectra from Gas Chromatographic Complex Mixture Analysis—Chemical Deconvolution

Kevin Siek<sup>1</sup>; Vasily Makarov<sup>2</sup>; Viatcheslav Artaev<sup>1</sup>; Dmitry Mazur<sup>3</sup>; Albert T. Lebedev<sup>3</sup> | <sup>1</sup>LECO Corporation, Saint Joseph, MI USA; <sup>2</sup>Mass Spectrometry Consulting Ltd., Bar, Montenegro; <sup>3</sup>Moscow State University, Moscow, Russian Federation

#### Introduction

- Impure spectra contain signals from multiple analytes or signals unrelated to the main analyte, due to propinquitous coelution that confounds purely mathematical de-coelution. Some such spectra are expected in complex sample analysis.
- Impure spectra may impede manual or automated interpretation and beget false assignments or other failures.
- information within high-resolution spectra can assist Chemical mathematical de-coelution and facilitate analyte identification.

# Methods

- Samples were prepared using conventional practices and analyzed by a Pegasus<sup>®</sup> GC-HRT (see handout for details).
- Lake Ontario fish tissue
- Municipal wastewater
- Nigerian crude oil
- The wastewater sample was subsequently analyzed by Pegasus GC-HRT 4D to confirm results (details in handout).
- Chromatographic feature finding was performed using commercially available software (ChromaTOF-HRT<sup>®</sup> 1.81).
- Peak tables were sent to chemical deconvolution software.

# Chemical Deconvolution (CD) Algorithm Steps

- Deisotope:
- Assign 2+ and 1+ charge states
- Assign provisional formulas to monoisotopic signals
- Assemble fragmentation graphs, constraining connections by:
- Chemically possible mass defect changes
- Rational RDBE, halogen count, and heteroatom count limits
- Other ad-hoc rules, such as forbid CH<sub>2</sub> loss
- Analyze each fragmentation graph in the spectrum separately:
- Remove "orphans" (one-point, no-edge fragmentation graphs)
- Find best molecular ion or top fragment formula for each graph
- Library search each graph and compare to full spectrum search
- Return results for each input spectrum:
- Elemental compositions of molecular ion(s) or top fragment(s)
- Merged library hits, filtered by minimum elemental composition

#### **Results Summary**

- False signals rejected from 40 hand-selected impure spectra
- 40 spectra with total of 64 spurious signals at high m/z and intensity
- 55 of these 64 spurious signals were orphaned
- 5 of these 64 were joined to a separate minor fragmentation graph
- 4 of these 64 joined to the principal analyte's fragmentation graph
- Fish tissue example coelution
- Incurred tetraCl-diMeO-benzene under fish fat and siloxane peaks
- CD separated analyte's fragmentation graph from matrix graphs
- Filtering library hits by formula yields only the one correct hit
- Wastewater example coelution
- Full spectrum reflects only benzothiazole, overlooking other analytes
- CD reports 3 analytes and 1 matrix component within this peak
- In GC-HRT 4D analysis, same 3 analytes separate in 2<sup>nd</sup> dimension
- Crude oil example coelution
- Coelution of 5 analytes sorted correctly
- CD generally finds chemically distinct analytes, as in the example
- CD generally fails for chemically similar analytes, as expected

# Results Example I: 3x Coelution in Fish Tissue Matrix

CD returns the one correct hit after filtering by formula; matrix (low m/z) + siloxane (high m/z) interference removed:





## **Results Example II: 4x Coelution in Wastewater**

**CD** distinguishes 3 molecular ions and 1 top fragment:



Same analytes found in subsequent GC-HRT 4D analysis:



### **Results Example III: 5x Coelution in Crude Oil**

Five separate molecular ions found in coelution at  $\approx$  969 s:



#### Conclusions

- Chemical deconvolution algorithms accurately reported independent components of dead coelutions where such components belong to distinct chemical classes.
- Present algorithms fail to distinguish chemically similar analytes such as branched and linear alkanes, thus chemical deconvolution is not a substitute for GCxGC
- Chemical deconvolution algorithms accurately rejected more than 80% of spurious signals from manually curated spectra found to be contaminated with unrelated signals.